



Philips Sonicare G3
Premium Gum Care
Standard sonic toothbrush
heads

3-pack

Standard size
Click-on
BrushSync mode pairing

HX9053/24

Healthier gums for a healthier smile

Up to 7x healthier gums in just 2 weeks*

Help patients practice better gum health and protect their smile with our G3 Premium Gum Care brush head. Soft sides flex and adapt to the contours of gums for 2x more surface contact* as they brush & superior plaque removal at the gumline

Innovative technology

- Philips Sonicare's advanced sonic technology
- Handle automatically selects optimal mode and intensity*

Designed to optimize performance

- Tested to meet your oral health needs
- For the most effective cleaning, every day
- Works with any Philips Sonicare click-on toothbrush

For exceptional gum care

- Up to 100% less gingival inflammation*

An unrivaled deep clean


- Up to 2x more surface contact* for an effortless deep clean
- Up to 10x more plaque removal than a manual toothbrush

PHILIPS

sonicare


Highlights

Focus on gingival health



With Premium Gum Care, even the deepest clean is gentle. As the brush head moves along the gumline, its flexible sides and bristles absorb any excessive pressure so gums are protected even if patients are brushing too hard. Its rounded profile allows the bristles to clean evenly while being gentle on gums. Extra-soft gumline bristles deliver a gentle yet effective clean along the gumline, where gum disease starts. The brush head's smaller size allows easy maneuvering even in hard-to-reach areas.


Adaptive cleaning technology



Patients get a personalized clean every time they brush with our adaptive cleaning technology. Soft, flexible rubber sides let Premium Gum Care adapt to the unique contours of the mouth. Our bristles adjust to gums and teeth, giving patients up to 2x more surface contact


than a DiamondClean brush head for deeper cleaning, even in hard-to-reach spots. Adaptive cleaning technology also: allows gentle tracking along the gumline; absorbs excessive brushing pressure; and enables an enhanced sweeping motion, for a unique in-mouth feel and superior cleaning.

Powerful plaque removal



Thanks to its flexible design, Premium Gum Care is proven to remove up to 10x more plaque than a manual toothbrush. It moves to enhance our unique Philips Sonicare cleaning technology so that however patients brush, they get an exceptional clean they can see and feel.

BrushSync mode pairing



Patients will always get their best possible clean and gum care with our BrushSync mode pairing feature. A microchip in the Premium Gum Care


brush head syncs with the connected Philips Sonicare toothbrush handle. It prompts the handle to automatically select the optimal brushing mode and intensity level for exceptional gum health. All your patients need to do is start brushing.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Smart brush head replacement



Brush heads exhibit fatigue and bristle wear after three months of normal use. Philips Sonicare smart toothbrushes provide accurate brush head replacement reminders based on how often patients brush and how much pressure they use. For those who don't use a smart toothbrush, Premium Gum Care's blue reminder bristles fade to white to let patients know it's time for a replacement brush head.

Standard sonic toothbrush heads
3-pack Standard size, Click-on, BrushSync mode pairing

HX9053/24

Specifications

Design and finishing

- Color: White

Design and finishing

- Smart Brush Head Recognition
- Bristle stiffness feel: Soft
- Material brush head: Soft, flexible rubber sides
- Reminder bristles: Blue bristle color fade away
- Size: Standard

Compatibility

- Brush head system: Click-on
- Suitable for these models: 2 Series plaque control, 2 Series plaque defense, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean,

Essence+, FlexCare, FlexCare Platinum, FlexCare Platinum Connected, FlexCare+, for Kids, HealthyWhite, HealthyWhite+, PowerUp

Items included

- Brush heads: 3 G3 Premium Gum Care

Quality and performance

- Replacement: Every 3 months
- Tested: for optimal usage

Health benefits

- Gum health: Up to 7x healthier gums*
- Plaque removal: Removes 10x more plaque****



Issue date 2022-05-23

Version: 2.2.1

EAN: 00 07502 00650 32

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* *in Gum Care Mode vs a manual toothbrush; measured by GBI
* *than a DiamondClean brush head
* **in Gum Care Mode vs a manual toothbrush; measured by GBI
* *** BrushSync™ mode-pairing is only compatible with Philips Sonicare
BrushSync™ enabled toothbrush handles
* ****than a manual toothbrush