

# PHILIPS

## AVENT

### Double electric breast pump

Comfort

Includes 2x 4oz Natural bottle



SCF334/25



## More milk and more comfort, in less time\*

### Breast pump with massage cushion

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

#### Easy to set up, express and clean

- Includes our Natural bottle and nipple for natural latch on
- Few separate parts and intuitive design

#### Effective expression with 3 customizable settings

- Choose the most effective and comfortable setting for you

#### More milk in less time\*

- Reduce expression time with double pumping

#### Soft massage cushion with massaging petals

- More comfortable pumping position due to unique design
- Soft massage cushion with massaging petals

# Highlights

## More milk in less time

Reduce expression time with double pumping compared to single pumping.

## Simple settings to choose from

When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

## Soft massage cushion

Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

## Relaxing position

The breast pump has a unique design, so your milk flows directly from your breast into the bottle or container, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to

lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

## Easy set up, express and clean

Few separate parts and intuitive design. Easily express milk for your baby, anytime. Closed system, no milk spillage in the tubes. Designed for hygienic expressing.

## Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

# Specifications

## Country of origin

England

## Development stages

Stage: 0 - 6 months

## Material

**Bottle:** BPA free\*, Polypropylene

**Nipple:** BPA free\*, Silicone

**Breast pump:** BPA free\* (food contact parts only)

## What is included

**Extra soft newborn flow nipple:** 2 pcs

**Sealing disc for milk storage:** 2 pcs

## Design

**Breast pump design:** Compact design

**Bottle design:** Ergonomic shape, Wide neck

## What is included

**Breast pad sample packs:** 2 pcs (2 sets of 2x Day pads & 2 sets of 2x Night pads) pcs

**Natural bottle 4oz:** 2 pcs

**Breast pump body:** 2 pcs

**Base unit incl. tubing:** 1 pcs

## Ease of use

**Breast pump use:** Easy cleaning, Fully compatible range, Intuitive assembly

## Functions

**No leaning forward:** Sit in a comfortable position

**Soft massage cushion:** Gentle stimulation

**Settings:** 3 Expression settings, 1 Stimulation mode

**More milk in less time:** The perfect time saver

## What is included

**Large massage cushion (25mm):** 2 pcs (A smaller size cushion is available separately)



\* More milk in less time: Reduce expression time with double pumping compared to single pumping.  
\* 60% of US breastfeeding mothers preferred the Philips Avent Comfort Double electric breast pump for comfort when expressing vs 29% of mothers who preferred the Medela in Style-Independent home placement test among 76 mothers in USA, August 2012. Medela is a registered trademark of Medela Holding AG  
\* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011