

# PHILIPS

## sonicare

2 Series plaque  
control

1 mode

1 brush head



HX6211/99



## Removes up to 7X more plaque\*

\*than a manual toothbrush

Offers you superior plaque removal at the push of a button

#### **Proven to improve oral health**

- Safe and gentle for your oral health needs

#### **Designed around you**

- Easy start program builds your Philips Sonicare routine
- Timers encourage thorough 2 minutes of brushing

#### **Provides a superior clean**

- Our unique technology brings you a powerful yet gentle clean
- Angled neck to help tackle cleaning in tough areas
- Remove up to 7x more plaque with our ProResults brush head
- Better plaque removal to help reduce cavities

# Highlights

## Remove up to 6x more plaque



Click on our ProResults gum health brush head for exceptional gum care. Our bristles give you an effective, gentle clean even in hard-to-reach spots. You get up to 7x better plaque removal\*, and gentle cushioning for an extra comfortable experience.

## Helps reduce cavities



Provides better plaque removal and access to back teeth; brushing twice a day every day with this brush can help reduce cavities

## Philips Sonicare technology



Powerful sonic vibrations whip up your toothpaste into plaque fighting bubbles, and drive them deep between your teeth and along


your gum line. At the same time your teeth will experience 62,000 gentle but effective brush movements. You'll get a month's worth of manual brushing in just 2 minutes.

## Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

## QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

## Reach those tricky spotsAngled



The unique shape of your ProResults gum health toothbrush handle combined with our brush heads means that even tricky spots like back teeth get a thorough clean.

## Safe and gentle



Whatever your specific needs, you can be sure of a safe brushing experience with the 2 Series plaque control. Our sonic technology is suitable for use with braces, fillings, crowns and veneers, and can even be used in the daily treatment of gum disease.

# Specifications

## Power

**Voltage:** 110-220 V

## Technical specifications

**Operating time (full to empty):** Up to 2 weeks

**Battery:** Rechargeable

**Battery type:** NiMH

## Design and finishing

**Color:** White on Guacamole

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Battery indicator:** Light shows battery status

**Brush head system:** Easy click-on brush heads

**Handle:** Slim ergonomic design

**Brushing time:** Up to 2 weeks

## Items included

**Charger:** 1

**Brush heads:** 1 ProResults plaque control

**Handle:** 1 2-series plaque control

## Cleaning performance

**Performance:** Removes up to 7x more plaque\*

**Health benefits:** Helps reduce cavities

**Speed:** Up to 31000 brush movement/min

**Timer:** Quadpacer and SmarTimer

## Modes

**Clean:** For exceptional everyday clean

