

PHILIPS

sonicare

Sonic electric
toothbrush

ProtectiveClean
4100

Built-in pressure sensor

1 cleaning mode

1 x BrushSync feature



HX6811/01



Cleaner teeth. Gently does it.

Remove up to 7x more plaque than a manual

Feel the difference of a gentle clean with our Pressure Sensor while removing up to 7x more plaque than a manual toothbrush

Proven to improve oral health

- Safe & gentle on sensitive areas, orthodontics & dental work
- Whitens teeth in just one week

Designed around you

- A great way to get used to your new toothbrush
- Encouragement to brush thoroughly

Innovative technology

- Always know when to replace your brush heads
- Philips Sonicare's advanced sonic technology
- Connects smart brush handle and smart brush heads
- Lets you know when you're pressing too hard

Highlights

Naturally whiter teeth



Click on the W2 Optimal White brush head to remove surface stains and reveal a whiter smile. With its densely packed central stainremoval bristles, it's clinically proven to whiten teeth in just one week.

Safe and gentle




You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns, and veneers, and helps prevent cavities and improve gum health.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Too much pressure?



You may not notice if you're brushing too hard, but your electric toothbrush with pressure sensor will. If you start taking it out on your gums, the toothbrush will make a pulsing sound as a reminder to ease off on the pressure.

BrushSync technology



A microchip-enabled technology that detects and synchronizes the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

BrushSync replacement reminder



All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to


replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

Easy-start program



Starting something new can take getting used to. That's why our easy-start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

Helpful timers



Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

Specifications

Power

Voltage: 110–220 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 2 weeks

Battery type: Lithium ION

Design and finishing

Color: Navy Blue

Service

Warranty: 2-year limited warranty

Ease of use

Battery indicator: Light shows battery status

Handle compatibility: Easy click-on brush heads

Handle: Slim ergonomic design

Brushing time: Up to 2 weeks

Items included

Charger: 1

Brush heads: 1 C2 Optimal Plaque Control

Handle: 1 ProtectiveClean

Cleaning performance

Performance: Removes up to 7x more plaque*

Health benefits: Helps reduce cavities

Speed: Up to 62000 brush movement/min

Timer: Quadpacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday clean

Smart sensor technology

BrushSync Replacement Reminder: Always know when to replace brush heads

Pressure sensor: Alerts when brushing too hard



* Removes up to 7x more plaque vs. a manual toothbrush
* Pressure sensor protects teeth and gums from excess brushing pressure
* Always know when to replace your brush head for an effective clean with BrushSync replacement reminder
* Quadpacer and Smart timer helps ensure dental professional recommended 2 minute brushing time and encourages brushing in each quadrant of the mouth
* Battery charge indicator light let's you know when to recharge; 2 week battery life