

Philips Sonicare Power
Flosser 3000
Cordless water flosser
with accessories

Removes up to 99.9% of plaque

Complete clean in 60 seconds
2 flossing modes, 3 intensities
250 ml reservoir



HX3806/33

Flossing reinvented with QuadStream technology

Removes up to 99.9% of plaque* in treated areas

Effortlessly thorough cleaning between teeth to improve your gum health. Unique Quad Stream tip cleans more surface area with less effort for more effective flossing every time, so it's easy to get flossing right.

Easy and convenient

- Easy-to-fill 250 ml reservoir
- Effortless technique
- Standard tip for spot cleaning
- Fast, universal charging with USB-C cable

Faster, more effective flossing

- 2 flossing modes, 3 intensities
- Complete clean in 60 seconds
- Pulse Wave technology guides from tooth to tooth
- Quad Stream technology for faster, more effective flossing

PHILIPS
sonicare

Cordless water flosser with accessories

Removes up to 99.9% of plaque Complete clean in 60 seconds, 2 flossing modes, 3 intensities, 250 ml reservoir


Highlights

Quad Stream technology



Unique X-shaped Quad Stream tip separates flow into 4 water jets, which cover more surface area between teeth and along the gum line. Get a faster, deeper clean without the hassle of string floss.

Pulse Wave technology



Gentle pulses of water keep you going and guide you from tooth to tooth in Deep Clean mode, so you get it right every time.

2 modes, 3 intensities



Choose the clean that suits your smile. Clean mode delivers continuous flow for effective everyday cleaning. Deep Clean pulses for more

thorough cleaning. Adjustable intensity level adapts for comfort.

Complete clean in 60 seconds



Access those hard-to-reach areas with a nozzle that rotates 360 degrees, so that you can use the device in any orientation. Get a thorough 360° clean in just 60 seconds from start to finish.

USB-C Charging



Convenient charging with a USB-C cable uses the same connection as most mobile phones or tablets. Get up to 14 days of cleaning from a single charge.

250 ml reservoir




Reservoir capacity of 250 ml holds enough water for a recommended 60 second clean. Simply twist and remove or use the side-fill door for easy refilling.

Effortless technique



Just point, press and clean for more effective flossing. Aim the tip between your teeth and along the gumline, and let it do all the work for you!

Standard tip for spot cleaning



A standard single-stream tip boosts pressure for spot cleaning and removing food debris when you need it.

Cordless water flosser with accessories

Removes up to 99.9% of plaque Complete clean in 60 seconds, 2 flossing modes, 3 intensities, 250 ml reservoir

Specifications

Power

- Voltage: Multi-voltage charger

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

Performance

- Cleaning: Clean entire mouth in 60–90 sec

Items included

- Cordless Power Flosser: 1
- F1 Standard nozzle: 1

- F3 Quad Stream nozzle: 1
- Travel bag: 1
- USB charging cable: 1
- USB wall adapter: 1

Modes

- Clean: For exceptional everyday cleaning
- Deep Clean+: For an invigorating deep clean
- Intensities: 3

Health benefits

- Plaque removal: Removes up to 99% of plaque in treated areas*

Design and finishing

- Colour: Black



Issue date 2023-10-12

Version: 11.9.4

EAN: 87 10103 99144 1

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* In an in-vitro study. Actual results may vary.