



Philips goLITE BLU  
Energy light

### Intense blue light

Cordless, travel lock & pouch  
as effective as 10,000 lux  
Treatment time indicator



HF3431/60

## Feel more energized with blue light, naturally

### Compact blue light energizes like bright daylight

Philips goLITE BLU is a portable energy light that uses the natural power of daylight to improve energy levels and mood. Just 20-30 minutes a day is all it takes to feel more energized and active, helping you to perform whenever you need to

#### **Advanced Philips LEDs provide comfortable light**

- Has a similar effect on well-being as a sunny blue sky
- Even brightness, ultimate eye comfort
- 100+ years of Philips lighting and healthcare know-how
- 100% UV free light - safe for eyes and skin

#### **Clinically proven: energizes naturally with light**

- Improves energy levels and mood
- Independent research confirms improved energy levels
- Fights energy dips, fatigue and winter blues
- Light is a natural energizer that fits a healthy lifestyle

#### **Compact on-the-go design, use it anywhere, anytime**

- Rechargeable battery and mains cord for use anywhere
- Light and durable go-anywhere design with protective pouch
- Results in just 20-30 min/day; use while reading, working
- As effective as much larger 10,000 lux white lights

# PHILIPS

## Energy light

Intense blue light Cordless, travel lock & pouch, as effective as 10,000 lux, Treatment time indicator

# Highlights

### Increases energy levels

Today's busy modern lifestyles demand peak performance all day long – but even the fittest person suffers energy crashes or lows, particularly when deprived of natural daylight, such as in an office in winter. goLITE BLU fights energy dips and improves mood using the natural energizing power of sunlight. With just a single 20-30 minute session proven to be effective, use goLITE BLU whenever you need to charge up and energize: to get going in the morning, reboot after lunch, or maintain sustained energy levels throughout the day.

### Light is a natural energizer

Daylight is essential to well-being and keeping fit and energized throughout the day – just like healthy nutrition and regular exercise. Philips goLITE BLU mimics the natural energizing power of daylight on a bright sunny day, stimulating special receptors in the eye that trigger the body's natural response to sunlight. This helps you to feel more energetic and active. On days when your energy level is low and you're feeling tired, goLITE BLU promotes vitality, naturally.

### Just 20-30 minutes a day

Just 20-30 minutes a day with goLITE BLU is proven to revitalize and promote a positive mood. Simply position the energy light at arm's length so that the light reaches your eyes. There's no need to look directly into the light, so you can multitask with goLITE BLU, benefitting from it at breakfast, or when working or reading. Since the light is UV free and comfortable and safe for eyes, you can use it for as long as you wish during the day (if used in the evening or before bedtime it may delay sleep). With daily use, you are likely to feel goLITE BLU's energizing effect within 1-2 weeks.

### Rechargeable battery



Rechargeable battery and mains cord for use anywhere.

### Ultimate eye comfort



EnergyUp uses a unique lighting system to control glare for complete eye comfort. A unique combination of diffusers, reflectors and filters distribute the powerful LED light evenly across the whole surface of the screen. This controls glare and prevents bright spots, for diffuse light that is always comfortable and pleasant to experience.

### UV-free

The pure blue light of goLITE BLU increases energy levels naturally. The effect is similar to that of a sunny day, but without the UV rays associated with sunlight.

HF3431/60

### Developed by Philips



Philips was founded more than one hundred years ago as a pioneer of affordable light bulbs, and later used its expertise in vacuum tube technology to develop one of the world's first X-ray imaging machines. To this day, the company remains a global leader in lighting and healthcare, driven on by the desire to improve people's lives around the world every day by creating innovations that matter to you. This proud heritage continues with goLITE BLU, an innovative product which benefits from all of Philips' experience as a leading lighting, healthcare and consumer technology company.

### Compact portable design



goLITE BLU's blue light is the most efficient kind of energizing light, requiring just 200 lux to achieve an effect similar to 10,000 lux of white light. Blue light's high efficiency helps ensure that goLITE BLU is our smallest energizing light yet. With a lightweight design, rechargeable battery and sturdy travel pouch, you can take this energy light with you anywhere, while it fits any desk or table at home or in the office too. Simply place goLITE BLU at arm's length within your field of vision and adjust the stand so that the light reaches your eyes.

## Energy light

Intense blue light Cordless, travel lock & pouch, as effective as 10,000 lux, Treatment time indicator

HF3431/60

# Specifications

### Weight and dimensions

- Product dimensions: 14,3 x 14,3 x 3,5 cm
- Product weight: 0,4 kg
- F-box weight: 0,8 kg
- F-box dimensions: 18,0 x 18,9 x 10,4 (WxHxD) cm
- Master carton: 2

### Logistic data

- Country of origin: China
- CTV code: 884343160

### Easy to use

- On/off button
- Treatment time indicator
- Light intensity memory
- Touch button operation
- Light intensity settings: 5

### Comfortable light

- Even screen luminance
- Specific blue light
- Similar effect as 10,000 lux

### Portable

- Small, compact design
- Protective pouch included

- Corded and cordless use
- Fully recharged battery: Up to 60 min treatment time
- Travel lock

### Well-being by light

- Fight winter blues
- Increase energy level
- Keeps you sharp
- Lift your mood
- Feel more rested

### Technical

- Voltage: 100/240 V
- Frequency: 50/60 Hz
- Insulation Power Plug: Class II (double isolation)
- Type of Lamp: High power LED
- Cordlength: 180 cm
- Lifetime of lamps: 10,000 hours
- Power adapter: 10W

### Safety and Regulations

- UV-free: No UV or near UV radiation
- Complies to Portable: Electric Luminaire UL153, CSA22.2 no 12
- cULus certified



Issue date 2023-09-21

Version: 2.2.1

UPC: 0 75020 04506 5

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* References:

\* [1] W.B. Duijzer, Y. Meesters (SLTBR 2011, p.33)