

Philips Sonicare Elite
Sonic electric toothbrush

HX9882



Custom Care system

4 ways of customized cleaning in one product

No two people are the same - especially when it comes to oral care. So why not use a power toothbrush that can adapt to your different needs and still deliver superior results.

2 brush head sizes and 2 cleaning speeds

- Thorough everyday clean, massages and stimulates gums
- Standard for a complete clean, compact for precision clean

Superior cleaning technology

- Clinically proven safe and gentle
- Dynamic fluid action reaches deep between teeth
- Removes plaque in hard to reach areas
- Healthier gums

Customizable cleaning

- Easy-start to acclimate to the Sonicare experience
- Quadpacer 30-second interval timer encourages even brushing

PHILIPS
sonicare


Highlights

2 brush head sizes




Standard for a complete clean, compact for precision clean

2 speeds




Max for a thorough everyday clean, Gentle for sensitive areas and gum massage.

Quadpacer



Quadpacer 30-second interval timer encourages even brushing

Easy-start



Gently increases power over first 14 uses, to ease into the Sonicare experience.

2-minute timer




Dynamic fluid action reaches deep between teeth

Deluxe recharge gauge

Removes plaque in hard to reach areas

Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

Healthier gums

Clinically proven to reduce gingivitis and help shrink periodontal pockets.

Specifications

Accessories

- HX7001: Standard Elite single pack brush head
- HX7002: Standard Elite double pack brush head
- HX7011: Mini Elite single pack brush head
- HX7012: Mini Elite double pack brush head

Technical specifications

- Voltage: 220-240 V
- Charging time: Hours hour(s)
- Operating time (full to empty): 28 brushings
- Battery rechargeable: NiCd
- Speed (top): Up to 31,000 movements per minute

Weight and dimensions

- F-box dimensions: 17.2 x 9.5 x 23.5 in
- F-box weight: 0.68 lb
- A-box dimensions: 35.1 x 29.7 x 24.1 in
- A-box weight: 4.08 lb
- A-box number of pieces: 6 pcs

Additional features

- Smartimer: Dentists' recommend a 2 minute brushing time
- Deluxe recharge gauge: Green tracks battery charge, yellow for recharge
- Ergonomic tapered handle: Non-slip grip designed to improve handling
- Easy-Start: Gently increases power over 14 uses
- Programmable Quadpacer: 30-second interval timer for an even clean
- 2 speeds: Choose the speed most comfortable for you
- 2 contoured brush head sizes: Standard for complete clean, Compact for precision

Cleaning performance

- Speed: Up to 31000 brush movement/min



Issue date 2022-01-17

Version: 2.2.1

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com