

# PHILIPS

## sonicare

Sonic electric  
toothbrush

HealthyWhite

3 Modes

HX6780/02



## Whiter, healthier teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. HealthyWhite is proven to remove everyday stains in just two weeks, with the regular use of the Clean & White Mode.

### **Proven to improve oral health**

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology

### **Guides you to follow dental recommendations**

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

### **Provides a superior clean**

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

### **Offers a customized brushing experience**

- Easy-start program to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2 minute mode
- Sensitive Mode: Gentle teeth and gum cleaning

# Highlights

## Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

## Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

## SmarterTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

## Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.\*

## Clean Mode

For exceptional everyday clean. Helps maintain whitening treatments.

## Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums



# Specifications

## Design and finishing

**Color:** Frost White

## Power

**Voltage:** 110-220 V

## Technical specifications

**Operating time (full to empty):** Up to 3 weeks\*\*

**Battery:** Rechargeable

**Battery type:** Lithium ION

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Handle:** Slim ergonomic design, Rubber grip for easy handling

**Brush head system:** Easy click-on brush heads

**Battery indicator:** Light shows battery status

**Brushing time:** Up to 3 weeks\*\*

## Items included

**Charger:** 1

**Handles:** 1 HealthyWhite

**Brush heads:** 1 ProResults standard

## Cleaning performance

**Performance:** Removes up to 90% of stains\*

**Health benefits:** Helps improve gum health

**Speed:** Up to 62000 brush movement/min

**Whitening benefits:** Whitens teeth up to 2 shades

**Timer:** Quadpacer and SmarTimer

## Modes

**Clean:** For exceptional everyday clean

**Sensitive:** Gentle teeth and gum cleaning

**Clean & White:** 2 min Clean + 30 sec White

