



Philips
EnergyLight



HF3318

Feel energetic

Boosts your energy as on a sunny day

The Philips Energy Light is a special light therapy device that gives you the daylight and energy for you to retain the summer feeling all year long in a way that easily fits your daily routine.

Safe to use

- UV-free

Feel energetic in 30 minutes a day

- 10,000 lux light intensity
- Clinically proven

Comfortable and easy to use light

- Adjustable angle
- Light intensity dimmer
- Treatment time indicator



PHILIPS

Highlights

10,000 lux light intensity



The EnergyLight produces light with all the intensity and quality of natural daylight. In fact, it can give you up to 10,000 lux, but without any unpleasant dazzle or glare.

Clinically proven



The Philips EnergyLight is clinically proven to increase your energy level and help fight winter blues.

Treatment time indicator



The EnergyLight has an LED timer which indicates the time the product has been used.

Light intensity dimmer




Fight winter blues and get the energy as you do in summer while using the appliance at a convenient level of light you for example need for easy-reading. The EnergyLight features a light intensity dimmer so it can easily be adjusted between zero to about 10,000 lux (at 15-20 cm) matching the intensity of your personal preference.

UV-free



The EnergyLight uses the full natural daylight spectrum. It meets European safety standards and features a UV-filtering screen. The EnergyLight is therefore safe to use.

Adjustable angle



The appliance can be positioned at various angles for comfortable use.

Specifications

Technical specifications

- Cord length: 200 cm
- Power: 75 W
- Frequency: 60 Hz
- Insulation: Class II (double insulation)
- Voltage: 120 V
- Type of lamps: 2 PL-L 36 W Philips EnergyLight
- Lifetime of lamps: 1,700 hours (10,000 lux)
- Ballast operating frequency: 42 kHz...125 kHz

Weight and dimensions

- F-box dimensions (=A-box): 36x 62.2 x 21.5 (W x H x D) cm
- F-box weight (=A-box): 4.4 kg

Logistic data

- Country of origin: Germany
- CTV code: 884331860

Safety

- UV-free: No UV or near UV radiation
- CE 0344 Certified

Easy to use

- On/off button
- Dimmer
- Treatment time indicator

Comfortable light

- 10,000 lux
- Even screen luminance

Easy positionable

- Adjustable angle: 20 degrees

Well-being by light

- Fight winter blues
- Increase alertness
- Increase energy level



Issue date 2022-01-16

Version: 2.2.1

EAN: 00 07502 00144 29

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com