

PHILIPS

sonicare

Sonic electric
toothbrush

HealthyWhite

3 modes

2 brush heads

Whitening pen

HX6732/76



Whiter, healthier teeth

The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Sonicare HealthyWhite electric toothbrush is proven to remove everyday stains in just two weeks, with the regular use of the Clean & White Mode.

Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology
- Improves gum health in only two weeks

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- Removes up to 7x more plaque than a manual toothbrush

Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Sensitive Mode: Gentle teeth and gum cleaning


Highlights

Up to 7x better plaque removal




Removes up to 7x more plaque than a manual toothbrush for a superior clean

Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

Naturally whiter teeth




The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck




The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

Easy-start program




Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience



Highlights

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

SmarterTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time



Specifications

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 3 weeks**

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: Frost White

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design, Rubber grip for easy handling

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Brushing time: Up to 3 weeks**

Items included

Handles: 2 HealthyWhite

Brush heads: 2 ProResults plaque control

Travel case: 2

Charger: 2

Whitening pen: 2

Cleaning performance

Health benefits: Improves gum health in only two weeks

Performance: Removes up to 6X more plaque*, Removes up to 90% of stains*

Whitening benefits: Whitens teeth up to 2 shades

Timer: Quadpacer and SmarTimer

Speed: Up to 62000 brush movement/min

Modes

Clean: For exceptional everyday clean

Sensitive: Gentle teeth and gum cleaning

Clean & White: 2 min Clean + 30 sec White

