



Philips Sonicare  
TongueCare+  
Tongue brushes

**2 pack**

Click-on  
Latex free  
BrushSync mode pairing

**HX8072/01**

## A Philips Sonicare clean, for your tongue

The only tongue cleaner with Sonicare strength

Transform a Philips Sonicare toothbrush into a sonic-powered tongue cleaner. With specially designed MicroBristles, our smart TongueCare+ tongue brush deep cleans your tongue when paired with Philip Sonicare BreathRx tongue spray.

### **An easy addition to your oral care routine**

- From toothbrush to tongue cleaner in just one click
- Easily reach every part of your tongue

### **Innovative technology**

- 240 MicroBristles for a gentle deep clean
- A Philips Sonicare clean for your tongue
- Automatically pairs with your toothbrush handle

**PHILIPS**

**sonicare**

Tongue brushes  
2 pack Click-on, Latex free, BrushSync mode pairing


## Highlights

### Easy to click-on




TongueCare+ tongue brush quickly transforms your Philips Sonicare toothbrush into a sonic-powered tongue cleaner. Each tongue brush clicks onto your Philips Sonicare toothbrush handle, just like a regular brush head. Simple to use, replace and clean, these innovative tongue brushes are compatible with all Philips Sonicare click-on toothbrush handles, so it couldn't be easier to add tongue cleaning to your daily oral care routine. It fits all Philips Sonicare toothbrush handles except: PowerUp Battery and Essence.

### Get your tongue truly clean



To get your tongue truly clean, you need to be able to reach everywhere. Thanks to its ergonomic shape, the TongueCare+ tongue brush lets you do this in comfort. The small and flexible head makes it easy to clean away bad-breath bacteria in every nook and cranny.


### Smart mode-pairing



Our smart TongueCare+ tongue brush has a microchip inside that communicates with the


toothbrush handle. So when you click it on, the handle automatically recognizes it and selects the appropriate mode and intensity level to optimize your tongue-cleaning performance. All you have to do is to press the power button.

### MicroBristles



With 240 rubber MicroBristles, the TongueCare+ tongue brush is specifically designed to clean the soft and porous surface of your tongue. The flexible MicroBristles contour to and gently clean around the ridges and grooves, to remove bacteria build-up and drive our bacteria-killing Philips Sonicare BreathRx tongue spray deeper. The compact shape allows you to comfortably clean the entire tongue.

### Sonic technology



TongueCare+ is the only tongue cleaner with Philips Sonicare strength, to clean away bad-breath bacteria with every stroke. 31,000 vibrations per minute break up tongue debris, sweep away bad-breath bacteria, and drive the bacteria-killing ingredients of the Philips Sonicare BreathRx tongue spray deeper. Your tongue will feel cleaner and your breath will stay fresh longer.

HX8072/01

## Specifications

### Design and finishing

- Color: White

### Design and finishing

- Smart Brush Head Recognition

### Material

- Elastomer: Latex-free

### Compatibility

- Brush head system: Click-on
- Suitable for these models: 2 Series plaque control, 2 Series plaque defense, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean, Essence+, FlexCare, FlexCare Platinum, FlexCare Platinum Connected, FlexCare+, HealthyWhite, HealthyWhite+, PowerUp, Series 1100, 2100 Series, Series 3100, Prestige 9900, DiamondClean 9000, ProtectiveClean

### Items included

- Brush heads: 2 TongueCare+ tongue brushes



Issue date 2022-06-27

Version: 3.3.1

EAN: 00 07502 00654 76

© 2022 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)