




Philips Sonicare  
AirFloss Pro/Ultra -  
Interdental cleaner

w/6100 ProtectiveClean handle

1 nozzle  
1 brush head



HX8492/72

## Healthier gums in 2 weeks, guaranteed\*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra\*\*\*\* is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.\*\*

### Improves oral health

- Clinically proven as effective as floss for gum health\*\*
- Up to 99.9% plaque removal\*\*\*
- Helps prevent cavities between teeth

### Provides a superior clean

- Up to 7x more plaque removal than a manual toothbrush
- Three modes, three intensity settings

### Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

### The easiest way to effectively clean between teeth

- Total daily confidence of a fresh and clean mouth
- Takes just 60 seconds: Point, Press, Clean!
- An easy way to start a healthy habit

**PHILIPS**  
**sonicare**

# Highlights

## Say goodbye to plaque




Densely-packed, high-quality bristles remove up to 7x more plaque than a manual toothbrush.

## Brushing your way




This toothbrush allows you to customize your brushing with a choice of three modes and three intensities. Clean mode is standard for superior cleaning. White is the ideal mode for removing surface stains. And Gum Care mode adds an extra minute of reduced-power brushing, so you can gently massage your gums. Three intensities allow you to switch the order of 'higher' and 'lower' to align with the intensities from left to right.

## Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*\*

## Micro-droplet technology




Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

## Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss Pro/Ultra is an easy way to clean deeper between teeth, helping to form a healthy habit.

## Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.\*\* Helps improve gum health in as little as two weeks.

## Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

## Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

## High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

# Specifications

## Power

- Voltage: Multi-voltage charger

## Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

## Design and finishing

- Color: White with grey accents

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Battery indicator: Handle LEDs show charge level
- Handle: Slim ergonomic design
- Battery Life: 1-3 weeks\*\*\*\*
- Nozzle attachment: Easily snaps on and off

## Items included

- Brush heads: 1 W2 Optimal White
- Charger: 1
- Travel case: Soft travel case
- Handle: 1 ProtectiveClean
- AirFloss Pro/Ultra handle: 1
- AirFloss Pro/Ultra nozzle: 1

## Cleaning performance

- Health benefits: Cleans entire mouth in 60 sec.

Improve gum health in 2 weeks\*

- Performance: Removes up to 99.9% plaque\*\*
- Clinically proven results: As effective as floss, for gum health improvement, for plaque reduction, for reducing gum bleeding
- Customizable settings: Auto-burst or manual burst, Single/double/triple burst\*\*\*
- Whitening benefits: Whitens teeth 2X better\*
- For best results: Change nozzle every 6 months

## Modes

- 3 intensities: Low, Medium, High
- Clean: For exceptional everyday clean
- White: Removes surface stains
- Gum Care: Gently massages gums

## BrushSync Mode Pairing

- C2 Optimal Plaque Control: Pairs with Clean mode
- G2 Optimal Gum Care: Pairs with Gum Care mode
- W DiamondClean: Pairs with White mode

## Smart sensor technology

- BrushSync Replacement Reminder: Always know when to, replace brush heads
- BrushSync Technology: Connects smart handle and, smart brush head



Issue date 2024-04-21

Version: 1.1.1

UPC: 0 75020 08111 7

© 2024 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* or your money back

\*\* when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.

\*\*\* From the treated areas; In a lab study, actual in-mouth results may vary

\*\*\*\* AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.

\*\*\*\*\* depending on burst setting used