

# IHE Integration Statement forConnect



Vendor	Product Name	Version	Date
Forcare	forConnect	2021-1	2021-02-15

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client (ITI-1)	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
Cross-Enterprise User Assertion (XUA)	X-Service-User (ITI-40)	
Internet User Authorization (IUA)	Resource Server (ITI-72)	OAuth Bearer Token
Cross-Enterprise Document Sharing (XDS.b)	Document Source (ITI-41)	Document Addendum
		Document Replacement
	Document Consumer (ITI-18, ITI-43)	Document Transformation
		Document Metadata Update
	Document Administrator	On-Demand Documents
		Delayed Document Assembly
	Document Administrator	Update Metadata (ITI-57)
Cross-Enterprise Document Reliable Interchange (XDR)	Document Source (ITI-41)	Transmit Home Community ID
Cross-Enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Source (ITI-15)	Set of DICOM Instances
		PDF Report
	Imaging Document Consumer	Text Report
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Source (ITI-41, RAD-68)	WADO Retrieve (RAD-55)
		Set of DICOM Instances (RAD-55, RAD-69)
		PDF Report
		CDA Wrapped Text Report
	Imaging Document Consumer (RAD-55, RAD-69)	CDA Imaging Report with Structured Headings
Cross Community Access (XCA)	Responding Gateway (ITI-38, ITI-39)	
Patient Identifier Cross-referencing for MPI (PIX)	Patient Identifier Cross-Reference Consumer (ITI-9)	
Patient Identifier Cross-referencing HL7V3 (PIXV3)	Patient Identifier Cross-Reference Consumer (ITI-45)	
Basic Patient Privacy Consents (BPPC)	Content Creator	
Scanned Documents (XDS-SD)	Content Creator	
Cross-Enterprise Document Workflow Profile (XDW)	Content Creator	XDW Cross-Community with XDR
	Content Consumer	Document Import Option
	Content Updater	Document Import Option
		XDW Cross-Community with XDR
Cross Enterprise Basic eReferral Workflow Definition Profile (XBeR-WD)	Referral Requester	

Internet address for vendor's IHE information [www.forcare.com](http://www.forcare.com)

## Links to Standards Conformance Statements for the Implementation

HL7 [www.forcare.com](http://www.forcare.com)

**Links to general information on IHE**

In North America: [www.ihe.net](http://www.ihe.net)

In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)

In Japan: [www.ihe-j.org](http://www.ihe-j.org)

---

© 2006-2021 Forcare Holding B.V. All rights reserved.

---