

IHE Integration Statement forAudit



| Vendor | Product Name | Version | Date |
|---------|--------------|---------|------------|
| Forcare | forAudit | 2021-1 | 2021-02-15 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

| Integration Profile | Actor | Option |
|--|---|--------------------|
| Consistent Time (CT) | Time Client (ITI-1) | |
| Audit Trail and Node Authentication (ATNA) | Secure Application (ITI-19, ITI-1, ITI-20) | |
| | Audit Record Repository (ITI-20) | |
| Retrieve ATNA Audit Event (ATNA) | Audit Record Repository (ITI-81) | |
| Cross-Enterprise User Assertion (XUA) | X-Service-User (ITI-40) | Subject-Role |
| Internet User Authorization (IUA) | Authorization Client (ITI-71, ITI-72) | OAuth Bearer Token |
| | Authorization Server (ITI-71) | OAuth Bearer Token |
| | Resource Server (ITI-72) | OAuth Bearer Token |
| Cross-Enterprise Document Sharing (XDS.b) | Document Consumer (ITI-18, ITI-43) | |
| Enterprise User Authentication (EUA) | Client Authentication Agent (ITI-2, ITI-3, ITI-4) | |
| | Kerberized Server (ITI-4) | |

Internet address for vendor's IHE information www.forcare.com

Links to Standards Conformance Statements for the Implementation

| | |
|-------|--|
| HL7 | www.forcare.com |
| DICOM | www.forcare.com |

Links to general information on IHE

| | | |
|--|---|--|
| In North America: www.ihe.net | In Europe: www.ihe-europe.net | In Japan: www.ihe-j.org |
|--|---|--|

© 2006-2021 Forcare Holding B.V. All rights reserved.