

# IHE Integration Statement forAudit



Vendor	Product Name	Version	Date
Forcare	forAudit	2021-2	2021-05-17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client (ITI-1)	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
	Audit Record Repository (ITI-20)	
Retrieve ATNA Audit Event (ATNA)	Audit Record Repository (ITI-81)	
Cross-Enterprise User Assertion (XUA)	X-Service-User (ITI-40)	Subject-Role
Internet User Authorization (IUA)	Authorization Client (ITI-71, ITI-72)	OAuth Bearer Token
	Authorization Server (ITI-71)	OAuth Bearer Token
	Resource Server (ITI-72)	OAuth Bearer Token
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer (ITI-18, ITI-43)	
Enterprise User Authentication (EUA)	Client Authentication Agent (ITI-2, ITI-3, ITI-4)	
	Kerberized Server (ITI-4)	

Internet address for vendor's IHE information [www.forcare.com](http://www.forcare.com)

## Links to Standards Conformance Statements for the Implementation

HL7	<a href="http://www.forcare.com">www.forcare.com</a>
DICOM	<a href="http://www.forcare.com">www.forcare.com</a>

## Links to general information on IHE

In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>
--	---	--

© 2006-2021 Forcare Holding B.V. All rights reserved.