

# IHE Integration Statement forIndex



| Vendor  | Product Name | Version | Date       |
|---------|--------------|---------|------------|
| Forcare | forIndex     | 2020-3  | 2020-08-18 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

| Integration Profile                                 | Actor   | Option                                       |
|---|---|--|
| Consistent Time (CT)                                | Time Client (ITI-1)                                       |  |
| Audit Trail and Node Authentication (ATNA)          | Secure Application (ITI-19, ITI-1, ITI-20)                |  |
| Cross-Enterprise User Assertion (XUA)               | X-Service-Provider (ITI-40)                               | Subject-Role<br>PurposeOfUse                 |
|   | X-Service-User (ITI-40)                                   | Subject-Role<br>PurposeOfUse                 |
| Internet User Authorization (IUA)                   | Resource Server (ITI-72)                                  | OAuth Bearer Token                           |
| Cross-Enterprise Document Sharing (XDS.b)           | Document Registry (ITI-18, ITI-42)                        | Patient Identity Feed (ITI-8)                |
|   |   | Document Metadata Update (ITI-57)            |
|   |   | On-Demand Documents (ITI-61)                 |
|   |   | Delayed Document Assembly Reference ID       |
| Remove Metadata and Documents (RMD)                 | Document Registry (ITI-18, ITI-62)                        |  |
| Multi-Patient Queries (MPQ)                         | Document Registry (ITI-51)                                |  |
| Patient Demographics Query (PDQ)                    | Patient Demographics Supplier (ITI-21)                    |  |
| Patient Demographics Query HL7 V3 (PDQV3)           | Patient Demographics Supplier (ITI-47)                    | Continuation Option                          |
| Patient Identifier Cross-referencing for MPI (PIX)  | Patient Identifier Cross-reference Manager (ITI-8, ITI-9) |  |
| Patient Identifier Cross-referencing HL7 V3 (PIXV3) | Patient Identifier Cross-reference Manager (ITI-45)       |  |
| Document Metadata Subscription (DSUB)               | Document Metadata Notification Broker (ITI-52, ITI-53)    | Document Metadata Publish Recipient (ITI-54) |
|   | Document Metadata Publisher (ITI-54)                      |  |

Internet address for vendor's IHE information [www.forcare.com](http://www.forcare.com)

## Links to Standards Conformance Statements for the Implementation

|       |  |
|-------|--|
| HL7   | <a href="http://www.forcare.com">www.forcare.com</a> |
| DICOM | <a href="http://www.forcare.com">www.forcare.com</a> |

## Links to general information on IHE

In North America: [www.ihe.net](http://www.ihe.net)      In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)      In Japan: [www.ihe-j.org](http://www.ihe-j.org)

© 2006-2020 Forcare Holding B.V. All rights reserved.