

IHE Integration Statement forStore



Vendor	Product Name	Version	Date
Forcare	forStore	2021-2	2021-05-17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client (ITI-1)	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
Cross-Enterprise User Assertion (XUA)	X-Service-Provider (ITI-40)	Subject-Role
		PurposeOfUse
	X-Service-User (ITI-40)	Subject-Role
		PurposeOfUse
Internet User Authorization (IUA)	Resource Server (ITI-72)	OAuth Bearer Token
Cross-Enterprise Document Sharing (XDS.b)	Document Repository (ITI-41, ITI-42, ITI-43)	
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Document Repository (RAD-68)	
Remove Metadata and Documents (RMD)	Document Repository (ITI-86)	
Internet address for vendor's IHE information www.forcare.com		

Links to Standards Conformance Statements for the Implementation

HL7	www.forcare.com
DICOM	www.forcare.com

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
----------------------------------------------------------------	-----------------------------------------------------------------------	------------------------------------------------------------

© 2006-2021 Forcare Holding B.V. All rights reserved.