

# IHE Integration Statement forBridge



Vendor	Product Name	Version	Date
Forcare	forBridge	2020-3	2020-08-18

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client (ITI-1)	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
Cross-Enterprise User Assertion (XUA)	X-Service-Provider (ITI-40)	Subject-Role PurposeOfUse
	X-Service-User (ITI-40)	Subject-Role PurposeOfUse
Internet User Authorization (IUA)	Resource Server (ITI-72)	OAuth Bearer Token
Patient Demographics Query HL7 V3 (PDQV3)	Patient Demographics Supplier (ITI-47)	Continuation Option
	Patient Demographics Consumer (ITI-47)	
Patient Identifier Cross-Referencing (PIX)	Patient Identifier Cross-Reference Consumer (ITI-9)	
Patient Identifier Cross-Referencing HL7 v3 (PIXV3)	Patient Identifier Cross-Reference Consumer (ITI-45)	
Cross Community Access (XCA)	Initiating Gateway (ITI-38, ITI-39)	XDS Affinity Domain Option (ITI-18, ITI-43) On-Demand Documents
	Responding Gateway (ITI-38, ITI-39)	On-Demand Documents
Cross Community Access for Imaging (XCA-I)	Initiating Imaging Gateway (RAD-69, RAD-75)	
	Responding Imaging Gateway (RAD-69, RAD-75)	
Cross Community Patient Discovery (XCPD)	Initiating Gateway (ITI-55)	
	Responding Gateway (ITI-55)	
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer (ITI-18, ITI-43)	Basic Patient Privacy Enforcement
	Document Source (ITI-41)	Document Addendum Document Replacement Document Transformation
Cross-Enterprise Document Reliable Interchange (XDR)	Document Recipient (ITI-41)	Basic Patient Privacy Enforcement Accepts Limited Metadata
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer (RAD-69)	
Cross-Community Document Reliable Interchange (XCDR)	Initiating Gateway (ITI-80)	Basic Patient Privacy Enforcement
	Responding Gateway (ITI-80)	Basic Patient Privacy Enforcement

Internet address for vendor's IHE information [www.forcare.com](http://www.forcare.com)

## Links to Standards Conformance Statements for the Implementation

HL7	<a href="http://www.forcare.com">www.forcare.com</a>
DICOM	<a href="http://www.forcare.com">www.forcare.com</a>

## Links to general information on IHE

In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>
--	---	--

© 2006-2020 Forcare Holding B.V. All rights reserved.