

IHE Integration Statement forIndex



| Vendor | Product Name | Version | Date |
|---------|--------------|---------|------------|
| Forcare | forIndex | 2021-2 | 2021-05-17 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

| Integration Profile | Actor | Option |
|-----------------------------------------------------|-----------------------------------------------------------|----------------------------------------------|
| Consistent Time (CT) | Time Client (ITI-1) | |
| Audit Trail and Node Authentication (ATNA) | Secure Application (ITI-19, ITI-1, ITI-20) | |
| Cross-Enterprise User Assertion (XUA) | X-Service-Provider (ITI-40) | Subject-Role PurposeOfUse |
| | X-Service-User (ITI-40) | Subject-Role PurposeOfUse |
| Internet User Authorization (IUA) | Resource Server (ITI-72) | OAuth Bearer Token |
| Cross-Enterprise Document Sharing (XDS.b) | Document Registry (ITI-18, ITI-42) | Patient Identity Feed (ITI-8) |
| | | Document Metadata Update (ITI-57) |
| | | On-Demand Documents (ITI-61) |
| | | Delayed Document Assembly Reference ID |
| Remove Metadata and Documents (RMD) | Document Registry (ITI-18, ITI-62) | |
| Multi-Patient Queries (MPQ) | Document Registry (ITI-51) | |
| Patient Demographics Query (PDQ) | Patient Demographics Supplier (ITI-21) | |
| Patient Demographics Query HL7 V3 (PDQV3) | Patient Demographics Supplier (ITI-47) | Continuation Option |
| Patient Identifier Cross-referencing for MPI (PIX) | Patient Identifier Cross-reference Manager (ITI-8, ITI-9) | |
| Patient Identifier Cross-referencing HL7 V3 (PIXV3) | Patient Identifier Cross-reference Manager (ITI-45) | |
| Document Metadata Subscription (DSUB) | Document Metadata Notification Broker (ITI-52, ITI-53) | Document Metadata Publish Recipient (ITI-54) |
| | Document Metadata Publisher (ITI-54) | |

Internet address for vendor's IHE information www.forcare.com

Links to Standards Conformance Statements for the Implementation

| | |
|-------|------------------------------------------------------|
| HL7 | www.forcare.com |
| DICOM | www.forcare.com |

Links to general information on IHE

| | | |
|----------------------------------------------------------------|-----------------------------------------------------------------------|------------------------------------------------------------|
| In North America: www.ihe.net | In Europe: www.ihe-europe.net | In Japan: www.ihe-j.org |
|----------------------------------------------------------------|-----------------------------------------------------------------------|------------------------------------------------------------|

© 2006-2021 Forcare Holding B.V. All rights reserved.