

# IHE Integration Statement forIndex



Vendor	Product Name	Version	Date
Forcare	forIndex	2021-1	2021-02-15

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client (ITI-1)	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
Cross-Enterprise User Assertion (XUA)	X-Service-Provider (ITI-40)	Subject-Role PurposeOfUse
	X-Service-User (ITI-40)	Subject-Role PurposeOfUse
Internet User Authorization (IUA)	Resource Server (ITI-72)	OAuth Bearer Token
Cross-Enterprise Document Sharing (XDS.b)	Document Registry (ITI-18, ITI-42)	Patient Identity Feed (ITI-8)
		Document Metadata Update (ITI-57)
		On-Demand Documents (ITI-61)
		Delayed Document Assembly Reference ID
Remove Metadata and Documents (RMD)	Document Registry (ITI-18, ITI-62)	
Multi-Patient Queries (MPQ)	Document Registry (ITI-51)	
Patient Demographics Query (PDQ)	Patient Demographics Supplier (ITI-21)	
Patient Demographics Query HL7 V3 (PDQV3)	Patient Demographics Supplier (ITI-47)	Continuation Option
Patient Identifier Cross-referencing for MPI (PIX)	Patient Identifier Cross-reference Manager (ITI-8, ITI-9)	
Patient Identifier Cross-referencing HL7 V3 (PIXV3)	Patient Identifier Cross-reference Manager (ITI-45)	
Document Metadata Subscription (DSUB)	Document Metadata Notification Broker (ITI-52, ITI-53)	Document Metadata Publish Recipient (ITI-54)
	Document Metadata Publisher (ITI-54)	

Internet address for vendor's IHE information [www.forcare.com](http://www.forcare.com)

## Links to Standards Conformance Statements for the Implementation

HL7	<a href="http://www.forcare.com">www.forcare.com</a>
DICOM	<a href="http://www.forcare.com">www.forcare.com</a>

## Links to general information on IHE

In North America: [www.ihe.net](http://www.ihe.net)      In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)      In Japan: [www.ihe-j.org](http://www.ihe-j.org)

© 2006-2021 Forcare Holding B.V. All rights reserved.