

# IHE Integration Statement forConnect



Vendor	Product Name	Version	Date
Forcare	forConnect	2020-1	2020-02-19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client (ITI-1)	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
Cross-Enterprise User Assertion (XUA)	X-Service-User (ITI-40)	
Cross-Enterprise Document Sharing (XDS.b)	Document Source (ITI-41)	Document Addendum
		Document Replacement
	Document Consumer (ITI-18, ITI-43)	Document Transformation
		Document Metadata Update
Cross-Enterprise Document Reliable Interchange (XDR)	Document Source (ITI-41)	On-Demand Documents
		Delayed Document Assembly
Cross-Enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Source (ITI-15)	Update Metadata (ITI-57)
		Transmit Home Community ID
Cross-Enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Source (ITI-15)	Set of DICOM Instances
		PDF Report
		Text Report
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer	WADO Retrieve (RAD-55)
		Imaging Document Source (ITI-41, RAD-68)
	Imaging Document Consumer (RAD-55, RAD-69)	
Cross Community Access (XCA)	Responding Gateway (ITI-38, ITI-39)	CDA Imaging Report with Structured Headings
Patient Identifier Cross-referencing for MPI (PIX)	Patient Identifier Cross-Reference Consumer (ITI-9)	
Patient Identifier Cross-referencing HL7V3 (PIXV3)	Patient Identifier Cross-Reference Consumer (ITI-45)	
Basic Patient Privacy Consents (BPPC)	Content Creator	
Scanned Documents (XDS-SD)	Content Creator	
Cross-Enterprise Document Workflow Profile (XDW)	Content Creator	XDW Cross-Community with XDR
	Content Consumer	Document Import Option
	Content Updater	Document Import Option
		XDW Cross-Community with XDR
Cross Enterprise Basic eReferral Workflow Definition Profile (XBeR-WD)	Referral Requester	

Internet address for vendor's IHE information [www.forcare.com](http://www.forcare.com)

## Links to Standards Conformance Statements for the Implementation

HL7	<a href="http://www.forcare.com">www.forcare.com</a>
DICOM	<a href="http://www.forcare.com">www.forcare.com</a>

**Links to general information on IHE**

In North America: [www.ihe.net](http://www.ihe.net)

In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)

In Japan: [www.ihe-j.org](http://www.ihe-j.org)

© 2006-2020 Forcare Holding B.V. All rights reserved.

---