

PHILIPS

AVENT

Pump, store, feed &
Care all-in-one set

Express, store, feed and care

Includes soft massage cushion

2 Bottles and 10 storage cups

Nipple protectors and cream

SCD292/01



More comfort, more milk* Pump, store, feed & care Single electric breast pump all-in-one set

The Breastfeeding Starter set is ideal for moms who plan to express their breast milk. It has all you need to pump & store your breast milk, and feed your baby. It includes specialty products to help you care for your breasts and nipples.

Express

- More comfortable pumping position due to unique design
- Features a gentle stimulation mode and 3 pumping settings
- Soft massage cushion with warm feel

Feed

- Unique Petals for a soft, flexible nipple without collapse
- Natural latch on due to the wide breast shaped nipple

Store

- For use in fridge/freezer
- Easily pump in milk storage container

Care

- Protect sore nipples during breastfeeding
- Silky soft feel

Highlights

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

Soft massage cushion




Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let down.

Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

For use in fridge/freezer

The Philips Avent cups can be stored in the fridge or freezer and are dishwasher safe.

Easily pump in milk storage cup

With the milk storage cup adapter, you can easily pump in the milk storage container. No leakage, no spillage of your precious breast milk.

Protect sore nipples during breastfeeding

Philips Avent Nipple Protectors are only designed for use when you have sore or cracked nipples and should be used with health professional advice.

Silky soft feel



Silky soft topsheet and breathable, natural materials. Dermatologically tested.

Specifications

Design

Bottle design: Wide neck, Ergonomic shape

Breast pump design: Compact design

Country of origin

England

Material

Breast pump: BPA free*, Polypropylene

Bottle: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

Breast pads: Dermatologically tested, Silky soft topsheet, Natural breathable material

Storage cups: Polypropylene, BPA free*

What is included

Milk storage cup lid: 10 pcs

Milk storage cup adaptor: 2 pcs

Sealing disc for milk storage: 1 pcs

Natural bottle 4oz: 1 pcs

Extra soft Newborn Flow nipple: 1 pcs

Natural bottle 9oz: 1 pcs

Extra soft Slow Flow nipple: 1 pcs

Breast pump body: 1 pcs

Standard size cushion: 1 pcs

Base unit incl. tubing: 1 pcs

Milk storage cup (180ml/6oz): 10 pcs

Travel cover: 1 pcs

Breast pad sample packs: 1 pack (2pcs)

Standard Nipple Protector: 2 pcs

Bra clip: 1 pcs

Moisturizing Nipple Cream: 15 ml

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Settings: 1 Stimulation mode, 3 Expression settings

Breast care: Day, Night

Development stages

Stage: Pregnancy, 0-6 months

