



Philips Sonicare 3100 series
Sonic electric toothbrush

Sonic technology

Pressure sensor
QuadPacer and SmarTimer
Slim ergonomic design



HX3671/13

Goodbye manual toothbrush. Hello Sonic technology.

Up to 3x better plaque removal*

Sonic technology combined with our brush action gently removes plaque up to 3x better* than a manual toothbrush.

Designed around you

- Always know when to replace your brush heads
- 14-day long lasting battery life
- Ergonomic design makes the toothbrush easy to hold and use

Optimize your brushing

- Optimize your brushing with SmarTimer and QuadPacer

Easy transition to Sonicare

- Easy-Start to ease the transition

Protects your gums

- Built-in pressure sensor protects your teeth & gums

Provides superior clean

- Up to 3x better plaque removal vs a manual toothbrush*
- Our unique technology brings you a powerful yet gentle clean

PHILIPS
sonicare

Highlights

Advanced Sonicare Technology



Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two month's worth of manual brushing in just 2 minutes.** 31000 brush strokes per minute gently clean your teeth, break up plaque, and sweep it away for an exceptional daily clean.

Removes plaque 3x better*



The Sonicare electric toothbrush with advanced sonicare technology is clinically proven to remove plaque up to 3X better* than a manual toothbrush. It removes plaque from your teeth and along the gum line while protecting your gums.

Protects your gums



Built-in pressure sensor automatically detects pressure you apply, warns you, and reduces

the toothbrush vibrations automatically to help protect your gums. The toothbrush will make a pulsing sound as a reminder to ease off the pressure. 7 out of 10 people found this feature helped them become better brushers.

SmartTimer and QuadPacer



2 minute SmarTimer and 30-second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

Designed around you



Slim, light-weight, and ergonomic design makes the toothbrush easy to hold and use so that you can get an amazing clean without any strain.

BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

Easy-Start



Our Easy-Start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

Long lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

Specifications

Power

- Voltage: DC5V

Technical specifications

- Battery: Rechargeable
- Operating time (full to empty): 14 days****
- Battery type: Lithium ION

Design and finishing

- Color: White

Service

- Warranty: 2-year limited warranty

Ease of use

- Battery indicator: Light shows battery status
- Handle: Slim ergonomic design
- Brush head system: Easy click-on brush heads
- BrushSync Replacement Reminder: Replacement reminder, Always know when to, replace brush heads, Reminder icon lights up, to always ensure

best results

Items included

- Charger: 1
- Handle: 1 3100 Series
- Brush head: 1 C1 ProResults standard

Cleaning performance

- Performance: Removes up to 3x more plaque*
- Speed: Up to 31000 brush movement/min
- Timer: Quadpacer and SmarTimer
- Pressure feedback: Vibrates handle to alert user

Modes

- Clean: For exceptional everyday clean

Smart sensor technology

- BrushSync Replacement Reminder: Always know when to, replace brush heads
- Pressure sensor: Alerts when brushing too hard



Issue date 2023-09-17

Version: 11.11.1

EAN: 87 10103 99750 4

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* vs manual toothbrush for healthier teeth and gums

* Individual results may vary

** Data on file

**** based on two periods of two-minute brushings per day on standard mode