

# PHILIPS

## sonicare



Sonic electric  
toothbrush

3100 series

Sonic technology

Pressure sensor

QuadPacer and SmarTimer

Slim ergonomic design



HX3671/13

## Goodbye manual toothbrush. Hello Sonic technology.

Up to 3x better plaque removal\*

Sonic technology combined with our brush action gently removes plaque up to 3x better\* than a manual toothbrush.

### **Designed around you**

- Always know when to replace your brush heads
- 14-day long lasting battery life
- Ergonomic design makes the toothbrush easy to hold and use

### **Optimize your brushing**

- Optimize your brushing with SmarTimer and QuadPacer

### **Easy transition to Sonicare**

- Easy-Start to ease the transition

### **Protects your gums**

- Built-in pressure sensor protects your teeth & gums

### **Provides superior clean**

- Up to 3x better plaque removal vs a manual toothbrush\*
- Our unique technology brings you a powerful yet gentle clean

# Highlights

## Advanced Sonicare Technology




Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two month's worth of manual brushing in just 2 minutes.\*\* 31000 brush strokes per minute gently clean your teeth, break up plaque, and sweep it away for an exceptional daily clean.

## Removes plaque 3x better\*




The Sonicare electric toothbrush with advanced sonicare technology is clinically proven to remove plaque up to 3X better\* than a manual toothbrush. It removes plaque from your teeth and along the gum line while protecting your gums.

## Protects your gums



Built-in pressure sensor automatically detects the pressure you apply, warns you, and reduces the toothbrush vibrations automatically to help protect your gums. The toothbrush will make a pulsing sound as a reminder to ease off the pressure. 7 out of 10 people found this feature helped them become better brushers.

## SmartTimer and QuadPacer



2 minute SmartTimer and 30-second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

## Designed around you



Slim, light-weight, and ergonomic design makes the toothbrush easy to hold and use so that you can get an amazing clean without any strain.

## BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

## Easy-Start



Our Easy-Start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

## Long lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

## 3100 series

Sonic electric toothbrush

HX3671/13

# Specifications

### Power

Voltage: DC5V

### Technical specifications

Battery: Rechargeable

Operating time (full to empty): 14 days\*\*\*\*\*

Energy Consumption: Standby without display  
<0.5 W

Battery type: Lithium ION

### Design and finishing

Color: White

### Service

Warranty: 2-year limited warranty

### Ease of use

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

BrushSync Replacement Reminder: Replacement reminder, Always know when to, replace brush heads, Reminder icon lights up, to always ensure best results

### Items included

Handle: 1 3100 Series

Brush head: 1 C1 ProResults standard

Charger: 1

### Cleaning performance

Speed: Up to 31000 brush movement/min

Performance: Removes up to 3x more plaque\*

Timer: Quadpacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

### Modes

Clean: For exceptional everyday clean

### Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

