

# PHILIPS

## AVENT

### Single electric breast pump

Pump anytime, anywhere

Includes soft massage cushion

Natural bottle and nipple



SCF332/31



## More comfort, more milk\*. Anytime, anywhere

### Single electric breast pump with massage cushion

Sit comfortably without leaning forward as our soft massage cushion gently stimulates your milk flow. Use our quiet pump anytime, anywhere, even with batteries. It's easy to set up, personalize, use and clean

#### Gently stimulates let down and milk flow

- Soft massage cushion with massaging petals

#### Easy to set up, express and clean

- Few separate parts and intuitive design
- Closed system, designed for hygienic expressing
- Includes our Natural bottle and nipple for natural latch on

#### Quiet pump for discreet expression, anytime

- Designed to let you express quietly

#### Clinically proven less leaning forward\*

- More comfortable pumping position due to unique design

#### 1 stimulation mode and 3 expression settings

- Choose the setting that is effective & comfortable for you

# Highlights

## Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle or container, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

## Simple settings to choose from

When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow comfortable for you.

## Soft massage cushion



Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

## No milk spillage in the tubes



No milk spillage in the tubes.

## Easy set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

## Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Quiet pump

Designed for discreet expression.

# Specifications

## Country of origin

England

## Development stages

Stage: 0 - 6 months

## Material

**Bottle:** BPA free\*, Polypropylene

**Nipple:** BPA free\*, Silicone

**Breast pump:** BPA free\* (food contact parts only)

## What is included

Sealing disc for milk storage: 1 pcs

## Design

**Breast pump design:** Compact design

**Bottle design:** Ergonomic shape, Wide neck

## What is included

Natural bottle 4oz: 1 pcs

Breast pump body: 1 pcs

Base unit incl. tubing: 1 pcs

Travel cover: 1 pcs

Breast pad sample packs: 1 pack (2pcs)

## Ease of use

**Breast pump use:** Easy cleaning, Fully compatible range, Intuitive assembly

## Functions

**No leaning forward:** Sit in a comfortable position

**Soft massage cushion:** Gentle stimulation

**Settings:** 1 Stimulation mode, 3 Expression settings

## What is included

**Standard size cushion(19.5 mm):** 1 pcs (a larger size cushion is available and sold separately)

**Disposable breast pads:** 4 pcs

