

# PHILIPS

## sonicare

Sonic electric toothbrush  
with case - pink

3100 series

Pressure sensor

2 minute SmartTimer

Up to 3 x better plaque removal

14 day long-lasting battery life



HX3673/11



## Goodbye manual toothbrush. Hello Sonic technology.

### Up to 3 x better plaque removal\*

Sonic technology combined with our brush action gently removes plaque up to 3 x better\* than a manual toothbrush.

#### **Designed around you**

- Always know when to replace your brush heads
- Makes traveling easier for you
- 14 day long-lasting battery life

#### **Optimise your brushing**

- Optimise your brushing with SmarTimer and QuadPacer

#### **Easy transition to Sonicare**

- Easy-Start to ease the transition

#### **Protects your gums**

- Built-in pressure sensor protects your teeth and gums

#### **Provides superior clean**

- Up to 3 x better plaque removal vs. a manual toothbrush\*
- Our unique technology brings you a powerful yet gentle clean


# Highlights

## Advanced Sonicare Technology



Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two months' worth of manual brushing in just 2 minutes.\*\* 31,000 brush strokes per minute gently clean your teeth, break up plaque and sweep it away for an exceptional daily clean.

## Removes plaque 3 x better\*



The Sonicare electric toothbrush with advanced sonicare technology is clinically proven to remove plaque up to 3 X better\* than a manual toothbrush. It removes plaque from your teeth and along the gum line while protecting your gums.

## Protects your gums



Built-in pressure sensor automatically detects the pressure you apply, warns you and reduces the toothbrush vibrations automatically to help

protect your gums. The toothbrush will make a pulsing sound as a reminder to ease off the pressure. 7 out of 10 people found this feature helped them become better brushers.

## SmartTimer and QuadPacer



2 minute SmartTimer and 30 second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

## BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep remind you when it's time to replace it.

## Easy-Start



Our Easy-Start programme gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

## Long-lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

## Slim travel case



Our premium travel case lets you store your toothbrush hygienically, while our compact USB charging base keeps you topped up when you are on the go. You get two weeks of regular use from a single full charge, but the charger's just the ticket for longer trips.

# Specifications

## Power

Voltage: DC5V

## Technical specifications

Operating time (full to empty): 14 days\*\*\*\*\*

Battery: Rechargeable

Battery type: Lithium ION

## Design and finishing

Colour: Sugar rose

## Service

Warranty: 2-year limited warranty

## Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

BrushSync Replacement Reminder:

Replacement reminder. Always know when to, replace brush heads, Reminder icon lights up, to always ensure best results

## Items included

Charger: 1

Travel case: 1

Handle: 1 3100 Series

Brush head: 1 C1 ProResults standard

## Cleaning performance

Speed: Up to 31,000 brush movements/min

Performance: Removes up to 3 x more plaque\*

Timer: Quadpacer and Smartimer

Pressure feedback: Vibrates handle to alert user

## Modes

Clean: For exceptional everyday cleaning

## Smart sensor technology

BrushSync Replacement Reminder: Always know when to, replace brush heads

Pressure sensor: Alerts when brushing too hard

