

PHILIPS

AVENT



Easy Comfort Single Electric Breast Pump

Includes one 120 ml/4 oz bottle

 Product Picture

SCF301/03

Convenience and Comfort for effective expression

Choose between 4 stimulation modes and 4 pumping settings to make every expression session comfortable and personalised. The lightweight, compact design also means expressing on-the-go is a breeze.

Lightweight and compact design

- Easy expressing, at home or on the go
- Battery and electricity operated

Multiple hassle-free settings for personalised use

- Hassle-free, set and forget settings
- 4 stimulation modes x 4 pumping modes for optimal expressing

Discreet and hygienic expressing

- Discreet sound and backlit display for discreet use, day/night
- Hygienic expressing

Rhythmic Wave Technology

- Designed to mimic babies' suckling rhythm

Designed to gently stimulate milk flow

- Soft massage cushion designed to help stimulate milk flow

Highlights

Personalised pumping

It's important to be comfortable when expressing and comfort means something different for every mum. With this breast pump, you can choose between 4 stimulation modes and 4 pumping settings. This way, you can make each expressing session efficient and comfortable.

Hassle-free settings

Mums have enough on their mind without having to recall their ideal setting - so we'll do that for you. Your preferred settings are stored, so at the press of a button you're ready to go. Need a break from pumping? Simply press the pause button, and to resume press play. It's that easy.

Soft massage cushion

Our breast pumps were designed with your comfort in mind. The soft and velvety texture of the massage cushions feel warm and comfortable against your skin, which can help stimulate milk flow. And the flexible material is designed to mimic your baby's suckling action, which can help milk letdown.

Lightweight and compact design

The lightweight and compact design means you can easily express at home or take it with you for expressing on the go.

Dual Power

The product runs on electricity. For extra mobility you can also insert batteries to allow cordless use. You can insert four non-chargeable 1.5 V AA batteries or four rechargeable 1.2 V AA batteries

Designed to mimic suckling

When babies breastfeed, they follow a rhythm of suck, swallow, breathe, with a ratio of 1:1:1. This breast pump features special Rhythmic Wave Technology that's designed to mimic a baby's suckling and swallowing pattern, pace and ratio (1:1), to quickly and efficiently generate milk flow. This means your expressing sessions feel natural and comfortable.

Discreet expressing

Our handy backlit display means you can express discreetly at night, without having to turn on the lights. And the soft sound means less disturbance.

Hygienic expressing

We know that every drop counts, that's why we make sure all breastmilk flows straight from breast to container, without any backflow into the tubes. All parts that come into contact with your breastmilk are made using BPA-free materials.

Specifications

Design

Bottle design: Ergonomic shape, Wide neck
Breast pump design: Compact design

Country of origin

China

Material

Bottle: BPA-free*, Polypropylene
Teat: BPA-free*, Silicone
Breast pump: BPA-free* (food contact parts only)

What is included

Sealing disc for milk storage: 1 pcs
Breast pump body: 1 pcs
Base unit incl. tubing: 1 pcs
Travel cover: 1 pcs
Breast pad sample packs: 1 pcs
Natural bottle 125 ml/4 oz: 1 pcs
Soft cushion: 1 pcs
Lanyard: 1 pcs

Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Functions

Rhythmic wave technology: Mimics baby suckling pattern

Development stages

Stage: 0–6 months

