



Philips Sonicare Elite
Sonic electric toothbrush

HX9882

Custom Care system

4 ways of customized cleaning in one product

No two people are the same - especially when it comes to oral care. So why not use a power toothbrush that can adapt to your different needs and still deliver superior results.

2 brush head sizes and 2 cleaning speeds

- Thorough everyday clean, massages and stimulates gums
- Standard for a complete clean, compact for precision clean

Superior cleaning technology

- Clinically proven safe and gentle
- Dynamic fluid action reaches deep between teeth
- Removes plaque in hard to reach areas
- Healthier gums

Customizable cleaning

- Easy-start to acclimate to the Sonicare experience
- Quadpacer 30-second interval timer encourages even brushing

PHILIPS
sonicare

Specifications

Accessories

- HX7001: Standard Elite single pack brush head
- HX7002: Standard Elite double pack brush head
- HX7011: Mini Elite single pack brush head
- HX7012: Mini Elite double pack brush head

Technical specifications

- Voltage: 220-240 V
- Operating time (full to empty): 28 brushings
- Charging time: Hours hour(s)
- Battery rechargeable: NiCd
- Speed (top): Up to 31,000 movements per minute

Weight and dimensions

- F-box dimensions: 17.1 x 9.50 x 23.5 in
- F-box weight: 0.78 lb
- A-box dimensions: 35.1 x 30.2 x 25.4 in
- A-box weight: 5.11 lb
- A-box number of pieces: 6 pcs pcs

Additional features

- Smartimer: Dentists' recommend a 2 minute brushing time
- Deluxe recharge gauge: Green tracks battery charge, yellow for recharge
- Ergonomic tapered handle: Non-slip grip designed to improve handling
- Easy-Start: Gently increases power over 14 uses
- Programmable Quadpacer: 30-second interval timer for an even clean
- 2 speeds: Choose the speed most comfortable for you
- 2 contoured brush head sizes: Standard for complete clean, Compact for precision

Cleaning performance

- Speed: Up to 31000 brush movement/min

Highlights

2 brush head sizes

Standard for a complete clean, compact for precision clean

2 speeds

Max for a thorough everyday clean, Gentle for sensitive areas and gum massage.

2-minute timer

Dynamic fluid action reaches deep between teeth

Deluxe recharge gauge

Removes plaque in hard to reach areas

Easy-start

Gently increases power over first 14 uses, to ease into the Sonicare experience.

Healthier gums

Clinically proven to reduce gingivitis and help shrink periodontal pockets.

Quadpacer

Quadpacer 30-second interval timer encourages even brushing

Safe and gentle

Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.



Issue date 2022-03-29

Version: 2.2.1

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com