

PHILIPS

 Product Picture

AVENT

Large massage
cushion

Fits all comfort breast pump

SCF167/01

Right size, maximum comfort

Breast pump with wider funnel

This soft massage cushion has a wider funnel, specially designed for mums with larger nipples. It can be used with Philips Avent Comfort breast pump range. See below images to help you decide if you need the larger breast pump cushion.

Comfortable pumping

- Wider funnel for comfortable fit
- Fits all Natural range breast pumps

Check your fit for optimal comfort

- Right fit: The nipple sits comfortably in the cushion funnel
- Choose large size: Nipple touches the walls of the funnel

Highlights

Large breast pump cushion

You can find out whether you need a large breast pump cushion by checking whether your nipple fits comfortably in the cushion funnel. If your nipple touches or rubs against the walls of the cushion funnel, then you should choose the large breast pump cushion to ensure a more comfortable expressing experience.

Correct fit

If the nipple sits comfortably in the cushion funnel, you are using the correct size massage cushion.

Too small fit

If your nipple touches or rubs against the wall of the cushion funnel, then you should choose the large breast pump cushion with wider funnel to ensure a more comfortable pumping experience.

Fits all Natural breast pumps

The large size soft massage cushion is compatible with all Philips Avent Comfort breast pumps in the Natural range.

Specifications

Country of origin

England

Development stages

Stage: Pregnancy, 0 - 6 months

Material

BPA free*

Replaceable part

Fits product types :: SCF330, SCF332, SCF334, SCF330, SCF332

What is included

Large size cushion: 1 pcs

Functions

Soft massage cushion: Gentle stimulation

Cushion

Material: Silicone, BPA free*

Compatibility

Compatible with Natural range

Ease of use

Cushion: Easy to clean, Dishwasher proof

