

PHILIPS

 Product Picture

AVENT

Comfort Manual
breast pump

Natural

Includes 4oz bottle

SCF330/11

More comfort, more milk naturally

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

Easy to combine breast and bottle feeding

- Includes our Natural bottle and nipple for natural latch on

Gently stimulates let down and milk flow

- Soft massage cushion with warm feel

Other benefits

- Intuitive assembly. Easy visual matching of parts
- Easy cleaning thanks to the small number of separate parts
- Easy manual operation
- Compatible with other Philips Avent feeding products

Easy to use, store and transport

- Compact lightweight design

Highlights

More comfortable position

The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Soft massage cushion

Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let down.

Includes our Natural nipple

The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Compact design

The breast pump has a compact design, which makes it easy to position. Its ergonomically shaped handle offers a comfortable hold and

full control when expressing milk. Being small and lightweight means it is easy to store and transport, making pumping on-the-go more discreet.

Intuitive assembly

Easy visual matching of parts for intuitive assembly.

Fully compatible range

This breast pump can be used in combination with other feeding products in the Philips Avent range, including our Classic bottles and milk storage containers. Philips Avent also provides a range of breast care accessories to help you breastfeed for longer and enhance your comfort.

Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. All parts are dishwasher proof.

Easy manual operation

The manual breast pump is ideal if you express milk occasionally and value compactness. Easily operated with one hand.

Specifications

Development stages

Stage: 0 – 6 months

Material

Breast pump: Polypropylene, BPA free*

Bottle: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

What is included

Extra soft newborn flow nipple: 1 pcs

Sealing disc for milk storage: 1 pcs

Design

Breast pump design: Compact design

Bottle design: Ergonomic shape, Wide neck

What is included

Breast pad sample packs: 2 (2 Day pads & 2 Night pads) pcs

Breast pump body with handle: 1 pcs

Natural bottle 4oz: 1 pcs

Standard size cushion: 1 pcs

Travel cover: 1 pcs

Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Breast pump

Material: BPA free*

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Functions

Includes Natural bottle: Combine breast&bottle feeding

