

IHE Integration Statement forConnect



| Vendor | Product Name | Version | Date |
|---|--|---|------------|
| Forcare | forConnect | 2022-1 | 2022-02-15 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below | | | |
| Integration Profile | Actor | Option | |
| Consistent Time (CT) | Time Client (ITI-1) | | |
| Audit Trail and Node Authentication (ATNA) | Secure Application (ITI-19, ITI-1, ITI-20) | | |
| Cross-Enterprise User Assertion (XUA) | X-Service-User (ITI-40) | | |
| Internet User Authorization (IUA) | Resource Server (ITI-72) | SAML Token | |
| Cross-Enterprise Document Sharing (XDS.b) | Document Source (ITI-41) | Document Addendum | |
| | | Document Replacement | |
| | Document Consumer (ITI-18, ITI-43) | Document Transformation | |
| | | Document Metadata Update | |
| Cross-Enterprise Document Reliable Interchange (XDR) | Document Source (ITI-41) | On-Demand Documents | |
| | | Delayed Document Assembly | |
| Cross-Enterprise Document Sharing for Imaging (XDS-I) | Imaging Document Source (ITI-15) | Update Metadata (ITI-57) | |
| | | Transmit Home Community ID | |
| Cross-Enterprise Document Sharing for Imaging (XDS-I) | Imaging Document Source (ITI-15) | Set of DICOM Instances | |
| | | PDF Report | |
| Cross-Enterprise Document Sharing for Imaging (XDS-I.b) | Imaging Document Source (ITI-41, RAD-68) | Text Report | |
| | | Imaging Document Consumer (RAD-55) | |
| Cross-Enterprise Document Sharing for Imaging (XDS-I.b) | Imaging Document Source (ITI-41, RAD-68) | Set of DICOM Instances (RAD-55, RAD-69) | |
| | | PDF Report | |
| | Imaging Document Consumer (RAD-55, RAD-69) | CDA Wrapped Text Report | |
| | | CDA Imaging Report with Structured Headings | |
| Cross Community Access (XCA) | Responding Gateway (ITI-38, ITI-39) | | |
| Patient Identifier Cross-referencing for MPI (PIX) | Patient Identifier Cross-Reference Consumer (ITI-9) | | |
| Patient Identifier Cross-referencing HL7V3 (PIXV3) | Patient Identifier Cross-Reference Consumer (ITI-45) | | |
| Basic Patient Privacy Consents (BPPC) | Content Creator | | |
| Scanned Documents (XDS-SD) | Content Creator | | |
| Cross-Enterprise Document Workflow Profile (XDW) | Content Creator | XDW Cross-Community with XDR | |
| | Content Consumer | Document Import Option | |
| | Content Updater | Document Import Option | |
| Cross Enterprise Basic eReferral Workflow Definition Profile (XBeR-WD) | Referral Requester | XDW Cross-Community with XDR | |
| Internet address for vendor's IHE information www.forcare.com | | | |

Links to Standards Conformance Statements for the Implementation

HL7 www.forcare.com

DICOM www.forcare.com

Links to general information on IHE

In North America: www.ihe.net

In Europe: www.ihe-europe.net

In Japan: www.ihe-j.org

© 2006-2022 Forcare Holding B.V. All rights reserved.
