

# IHE Integration Statement forView



Vendor	Product Name	Version	Date
Forcare	forView	2022-1	2022-02-15
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below			
Integration Profile	Actor	Option	
Consistent Time (CT)	Time Client (ITI-1)		
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)		
Enterprise User Authentication (EUA)	Client Authentication Agent Kerberized Server		
Cross-Enterprise User Assertion (XUA)	X-Service-Provider (ITI-40)	Subject-Role	
		PurposeOfUse	
	X-Service-User (ITI-40)	Subject-Role	
		PurposeOfUse	
Internet User Authorization (IUA)	Authorization Client (ITI-71, ITI-72)	JWT Token	
		SAML Token	
	Resource Server (ITI-72)	JWT Token	
		SAML Token	
Cross-Enterprise Document Sharing (XDS.b)	Document Source (ITI-41)	Document Replacement	
		Document Addendum	
		Document Transformation	
	Document Consumer (ITI-18, ITI-43)	Basic Patient Privacy Enforcement	
		Basic Patient Privacy Proof	
		Document Metadata Update	
Remove Metadata and Documents (RMD)	Document Administrator (ITI-18, ITI-62)	On-Demand Documents	
		Delayed Document Assembly	
		Document Administrator (ITI-57)	
Cross-Enterprise Document Reliable Interchange (XDR)	Document Source (ITI-41)	Transmit Home Community ID	
Cross-Enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Consumer (RAD-55, RAD-69)		
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer (RAD-55, RAD-69)		
Multi-Patient Queries (MPQ)	Document Consumer (ITI-51)		
Cross-Enterprise Document Reliable Interchange (XDR)	Document Source (ITI-41)		
Cross-Community Access for Imaging (XCA-I)	Imaging Document Consumer (RAD-69)		
Patient Demographics Query (PDQ)	Patient Demographics Consumer (ITI-21)		
Patient Demographics Query HL7 V3 (PDQV3)	Patient Demographics Consumer (ITI-47)		
Patient Demographics Query for Mobile (PDQm)	Patient Demographics Supplier (ITI-78)		
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer (ITI-9)		
Patient Identifier Cross-referencing HL7 V3 (PIXV3)	Patient Identifier Cross-reference Consumer (ITI-45)		

Patient Identifier Cross-reference for Mobile (PIXm)	Patient Identifier Cross-reference Manager (ITI-83)	
Personnel White Pages (PWP)	Personnel White Pages Consumer (ITI-24)	
Document Metadata Subscription (DSUB)	Document Metadata Notification Recipient (ITI-53)	
	Document Metadata Subscriber (ITI-52)	
Portable Data for Imaging (PDI)	Portable Media Creator (RAD-47)	USB Media Option
Basic Patient Privacy Consents (BPPC)	Content Creator	Basic Patient Privacy Acknowledgement
	Content Consumer	Basic Patient Privacy Acknowledgement with Scanned Document
		Basic Patient Privacy Acknowledgement View
Cross-Enterprise Sharing of Medical Summaries (XDS-MS)	Content Creator	Referral Option
	Content Consumer	View Option
Scanned Documents (XDS-SD)	Content Consumer	View Option
Sharing Laboratory Reports (XD-LAB)	Content Consumer	View Option
Cardiac Imaging Report Content (CIRC)	Content Consumer	View Option
Antepartum Education (APE)	Content Consumer	View Option
Antepartum History and Physical (APHP)	Content Consumer	View Option
Antepartum Record Profile (APR)	Content Consumer	View Option
Antepartum Summary (APS)	Content Consumer	View Option
Antepartum Laboratory (APL)	Content Consumer	View Option
Composite Triage and Nursing Note (CTNN)	Content Consumer	View Option
Emergency Department Encounter Summary (EDES)	Content Consumer	View Option
Emergency Department Physician Note (EDPN)	Content Consumer	View Option
Early Hearing Care Plan (EHCP)	Content Consumer	View Option
Emergency Department Referral (EDR)	Content Consumer	View Option
Healthy Weight (HW)	Content Consumer	View Option
Labor and Delivery History and Physical (LDHP)	Content Consumer	View Option
Labor and Delivery Summary (LDS)	Content Consumer	View Option
Immunization Content (IC)	Content Consumer	View Option
Newborn Discharge Summary (NDS)	Content Consumer	View Option
Nursing Note (NN)	Content Consumer	View Option
Patient Plan of Care (PPOC)	Content Consumer	View Option
Postpartum Visit Summary (PPVS)	Content Consumer	View Option
Triage Note (TN)	Content Consumer	View Option
Exchange of Personal Health Record Content (XPHR)	Content Consumer	View Option
Community Dispense (DIS)	Content Consumer	View Option
Community Prescription (PRE)	Content Consumer	View Option
Community Pharmaceutical Advice (PADV)	Content Consumer	View Option
Vital Records Death Reporting (VRDR)	Content Consumer	View Option

Physician Reporting to Public Health Repository - Cancer Registry (PRPH-Ca)	Content Consumer	View Option
Cross-Enterprise Document Workflow Profile (XDW)	Content Creator	
	Content Consumer	View Option
	Content Updater	View Option
Cross Enterprise Tumor Board Workflow Definition (XTB-WD)	TBR Requestor	Workflow Input and Output
	TBR Scheduler	Workflow Input and Output
	TBR Preparer	Workflow Input and Output
	TBR Report Writer	Workflow Input and Output
	TBR Finalizer	Workflow Input and Output
Cross Enterprise Basic eReferral Workflow Definition (XBeR-WD)	Referral Requester	
	Referral Scheduler	
	Referral Performer	Process without Scheduling Phase Option
Mobile access to Health Documents (MHD)	Document Responder (ITI-66, ITI-67, ITI-68)	
	Document Recipient (ITI-65)	
Internet address for vendor's IHE information <a href="http://www.forcare.com">www.forcare.com</a>		

#### Links to Standards Conformance Statements for the Implementation

HL7	<a href="http://www.forcare.com">www.forcare.com</a>
DICOM	<a href="http://www.forcare.com">www.forcare.com</a>

#### Links to general information on IHE

In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>
© 2006-2022 Forcare Holding B.V. All rights reserved.		