

IHE Integration Statement forNotify



Vendor	Product Name	Version	Date
Forcare	forNotify	2022-1	2022-02-15

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
Cross-Enterprise User Assertion (XUA)	X-Service-Provider (ITI-40)	Subject-Role PurposeOfUse
	X-Service-User (ITI-40)	Subject-Role PurposeOfUse
Internet User Authorization (IUA)	Resource Server (ITI-72)	SAML Token
Document Metadata Subscription (DSUB)	Document Metadata Notification Recipient (ITI-53)	
	Document Metadata Notification Subscriber (ITI-52)	
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer (ITI-18, ITI-43)	
Cross-Enterprise Document Workflow (XDW)	Content Consumer	Document Import Option

Internet address for vendor's IHE information www.forcare.com

Links to Standards Conformance Statements for the Implementation

HL7	www.forcare.com
DICOM	www.forcare.com

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--

© 2006-2022 Forcare Holding B.V. All rights reserved.