



TT3000/01





# Prevent hair regrowth

### enjoy smooth skin everyday

The Philips Lumea IPL hair removal system for men uses gentle pulses of light to stimulate the hair root into a resting phase. Hair sheds naturally and regrowth is inhibited. Using it every two weeks enables you to enjoy smooth skin every day.

#### Effective hair regrowth control

· Body hair regrowth control

#### Intense Pulsed Light (IPL) technology

· Suitable for home use

#### **Gentle application**

· Hairs are eliminated before they appear

#### **Application areas**

• For chest, stomach, legs, underarms, back and pubic area

#### Suitable hair and skin types

 $\bullet\,$  NOT suitable for red, grey, light-blond hair and skin type VI

### Highlights

#### **IPL** technology



Philips Lumea works with an innovative light-based technology called "Intense Pulsed Light" (IPL) derived from technology used in professional beauty salons for hair removal. Philips Lumea now brings this innovative new technology, suitable for home use, to be used by you wherever and whenever you want.

#### Suitable hair and skin types



Philips Lumea works effectively on (naturally) dark-blond, brown and black hair. Philips Lumea is NOT effective for (naturally) red, grey and light blond hair. Philips Lumea is NOT suitable for very dark skin (skin type VI)

#### Smoother results than a blade



Hairs are eliminated before they appear

#### **Application body areas**



Philips Lumea for Men has been developed in close co-operation with leading dermatologists for a safe, effective and gentle application on arms, underarms, shoulders, stomach, pubic area (excluding the scrotum), legs, chest and back

#### **Body hair regrowth control**



Philips Lumea for Men applies gentle pulses of light to the hair root. Consequently, the hair sheds naturally and hair regrowth is inhibited. Repeating this treatment every two weeks for the underarms, pubic area, chest, stomach, legs and back will keep your skin smooth every day. For the most effective treatment, preshave the application areas with the included Bodygroom for as long as you still experience regrowth

#### Lamp does not need replacement

Philips Lumea comes equipped with a durable lamp with a lifetime of 80,000 flashes

## Specifications

#### **Features**

3cm² application window: For an effective

treatment of large areas

5 light energy settings: Adjustable to your

personal skin type

Integrated UV filter: Protects your skin from

UV light

Ready-to-flash light: Indicates correct

application on the skin

#### **Accessories**

Storage pouch: included Cleaning cloth: included

AC/DC Adapter/Charger: included

#### **Technical specifications**

Battery type: Lithium Ion Frequency: 50–60 Hz

Power: 7.5 W Voltage: 100-240 V

Light energy level: 2-6.5 J/cm² (depending on

chosen setting)

Light spectrum: >570 nm

High performance lamp: Delivers > 80,000

flashes

#### Logistic data

Country of origin: Hungary

#### Items included

Instructions for use: Full user manual, Quick

Start Guide

#### Power

Charging time: 100 min

Power system: Rechargeable battery Running time: Min. 320 flashes

Voltage: 100-240 V

#### Safety and adjustable settings

Safety ring system: Prevent unintentional

flashing

#### Service

Warranty: 2 years global warranty

© 2019 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2019-07-30 Version: 6.1.1

