

# PHILIPS

## Microdermabrasion

### VisaCare

Skin resurfacing treatment  
with Sensitive and Normal tip



SC6240/01



## Firmer feeling, more youthful looking skin with Philips VisaCare Microdermabrasion System

Invigorate and maintain your skin's natural youthful look with VisaCare's Dual Action Air Lift & Exfoliation system. Massages the skin, encouraging blood circulation, while removing dead skin cells, allowing for natural cell renewal.

### **Firmer feeling, more youthful looking skin**

- For a more youthful looking skin
- Skin is more radiant, feels smoother and firmer in 4-6 weeks

### **Microdermabrasion technology**

- DualAction Air Lift & Exfoliation System
- The Air Lift System gently stimulates microcirculation
- The exfoliating tip stimulates the cell renewal process

### **Safe and gentle, clinically tested**

- Suitable for all skin types
- Developed with expert scientists and dermatologists
- Salon technology, for safe use in your own home

### **Easy to use**

- 5 minute treatment only, twice a week
- Cordless usage
- Charge overnight for 3 weeks of use




[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Youthful

Philips VisaCare Microdermabrasion stimulates the resurfacing process of the skin for smoother, firmer feeling, visibly radiant skin.

## Visible results



93% of women tested reported instantly smoother feeling skin. 84% agree their skin looks more radiant, 84% noticed their skin looks healthier\*.

## DualAction System




Massages the skin, encouraging blood circulation, while removing dead skin cells, allowing for natural cell renewal.

## Air Lift System

The vacuum massage intensively stimulates microcirculation. Skin is gently lifted and stretched, massaging its deeper layers. Blood flow increases, bringing oxygen and nutrients to the surface. For a revitalized look. VisaCare is designed to be in the sweet spot of stretch required to stimulate collagen and elastin production. For a firmer feel in 4-6 weeks time.

## Exfoliating tip



The exfoliating tip gently removes part of the stratum corneum (the upper layer of the skin) and thickens the epidermis (the second layer of the skin) unveiling skin that feels smoother and looks radiant. It gently nudges away roughness and flakiness more effectively and controlled than manually possible. Skin is stimulated to renew cells. For a firmer feel in 4-6 weeks time.

## Suitable for all skin types




VisaCare's gentleness is clinically proven. The device includes two tips (sensitive and normal) addressing different levels of skin sensitivity. Philips skincare experts are at your disposal to guide you in achieving the best results.

## Developed with skin experts


Philips VisaCare uses microdermabrasion technology commonly used in professional beauty salons. Philips has adapted this technology for safe use in the comfort of your own home. Philips VisaCare has been developed in collaboration with expert scientists and dermatologists.

## 5 minutes treatment



Achieve immediate results with only 5 minutes treatment twice a week. With a smooth and healthy glow!

## Cordless usage



VisaCare's unique cordless design allows you to enjoy your treatment in the most comfortable way.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Benefits

**Dead skin cells removal:** thanks to the exfoliating tip

**Epidermal thickening:** driven by abrasive system

**Stimulated microcirculation:** lead by Air Lift System™

**Skincare creams:** better absorbed, experienced by 87% of users

## Ease of use

**Customized tips:** Normal and sensitive skin

**Cordless:** up to 6 uses without charging

**Battery indicator:** illuminated icon

**Easy handling:** Ergonomic and elegant design

**Exclusive and convenient stand:** charging and storing stand

## Application areas

Face and neck

## Items included

**Customized tips:** Normal and sensitive skin

**Instructions for use:** Full user manual

**Exclusive and convenient stand:** charging and storing stand

**Power adapter:** 100 - 240 V adapter

**Storage pouch:** For storage and travel

**Cleaning brush:** To clean VisaCare tips

## Power

**Charging time:** 8 hours

**Power system:** Rechargeable battery

**Running time:** 6 uses of 5 minutes each

**Voltage:** 100-240 V

## Service

**Warranty:** 2-year limited warranty

