


PHILIPS

Microdermabrasion

VisaCare

Skin resurfacing treatment
with Sensitive and Normal tip

SC6240/00



Firmer feeling, more youthful looking skin

with Philips VisaCare Microdermabrasion System

Invigorate and maintain your skin's natural youthful look with VisaCare's Dual Action Air Lift & Exfoliation system. Massages the skin, encouraging blood circulation, while removing dead skin cells, allowing for natural cell renewal.

Firmer feeling, more youthful looking skin

- For a more youthful looking skin
- Skin is more radiant, feels smoother and firmer in 4-6 weeks

Microdermabrasion technology

- DualAction Air Lift & Exfoliation System
- The Air Lift System gently stimulates microcirculation
- The exfoliating tip stimulates the cell renewal process

Safe and gentle, clinically tested

- Suitable for all skin types
- Developed with expert scientists and dermatologists
- Salon technology, for safe use in your own home

Easy to use


- 5 minute treatment only, twice a week
- Cordless usage
- Charge overnight for 3 weeks of use

Highlights

Youthful

Philips VisaCare Microdermabrasion stimulates the resurfacing process of the skin for smoother, firmer feeling, visibly radiant skin.

Visible results



93% of women tested reported instantly smoother feeling skin. 84% agree their skin looks more radiant, 84% noticed their skin looks healthier*.

DualAction System




Massages the skin, encouraging blood circulation, while removing dead skin cells, allowing for natural cell renewal.

Air Lift System


The vacuum massage intensively stimulates microcirculation. Skin is gently lifted and stretched, massaging its deeper layers. Blood flow increases, bringing oxygen and nutrients to the surface. For a revitalized look. VisaCare is designed to be in the sweetspot of stretch required to stimulate collagen and elastin production. For a firmer feel in 4-6 weeks time.

Exfoliating tip



The exfoliating tip gently removes part of the stratum corneum (the upper layer of the skin) and thickens the epidermis (the second layer of the skin) unveiling skin that feels smoother and looks radiant. It gently nudges away roughness and flakiness more effectively and controlled than manually possible. Skin is stimulated to renew cells. For a firmer feel in 4-6 weeks time.

Suitable for all skin types




VisaCare's gentleness is clinically proven. The device includes two tips (sensitive and normal) addressing different levels of skin sensitivity. Philips skincare experts are at your disposal to guide you in achieving the best results.

Developed with skin experts


Philips VisaCare uses microdermabrasion technology commonly used in professional beauty salons. Philips has adapted this technology for safe use in the comfort of your own home. Philips VisaCare has been developed in collaboration with expert scientists and dermatologists.

5 minutes treatment



Achieve immediate results with only 5 minutes treatment twice a week. With a smooth and healthy glow!

Cordless usage



VisaCare's unique cordless design allows you to enjoy your treatment in the most comfortable way.

Overnight charge



Charge your VisaCare in the elegant charging stand and enjoy 3 weeks of treatment without the need to plug the device.

Salon technology

Salon technology, for safe use in your own home.

Specifications

Benefits

Dead skin cells removal: thanks to the exfoliating tip

Epidermal thickening: driven by abrasive system

Stimulated microcirculation: lead by Air Lift System™

Skincare creams: xx times better absorption

Ease of use

Customized tips: Normal and sensitive skin

Cordless: up to 6 uses without charging

Battery indicator: illuminated icon

Easy handling: Ergonomic and elegant design

Exclusive and convenient stand: charging and storing stand

Application areas

Face

Items included

Customized tips: Normal and sensitive skin

Instructions for use: Full user manual

Exclusive and convenient stand: charging and storing stand

Power adapter: 100 - 240 V adapter

Storage pouch: For storage and travel

Cleaning brush: To clean VisaCare tips

Power

Charging time: 12 hours

Power system: Rechargeable battery

Running time: 6 uses of 5 minutes each

Voltage: 100-240 V

Service

Warranty: 2-year limited warranty

