

Philips Sonicare  
AirFloss Pro/Ultra -  
Interdental cleaner

w/4500 ProtectiveClean handle

1 nozzle  
1 brush head



HX8424/31

## Healthier gums in 2 weeks, guaranteed\*

### Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra\*\*\*\* is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.\*\*

#### Improves oral health

- Up to 99.9% plaque removal\*\*\*

#### Optimize your brushing

- Improves gum health by up to 100%\*\*\*\*\*
- Lets you know if you're pressing too hard
- Easy start program builds your Philips Sonicare routine
- Always know when to replace your brush heads
- Philips Sonicare's advanced sonic technology

#### Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

#### The easiest way to effectively clean between teeth

- Total daily confidence of a fresh and clean mouth
- Takes just 60 seconds: Point, Press, Clean!

**PHILIPS**  
**sonicare**

# Highlights

## Improves gum health



Densely-packed, high-quality bristles give you an extra gentle brushing to remove plaque along the gumline to improve gum health up to 100% more than a manual toothbrush. And, the specially curved power tip makes reaching the teeth at the back of your mouth a breeze.

## Monitor your brushing pressure




Brushing too hard can damage your teeth and gums. To prevent this, your Philips Sonicare ProtectiveClean emits a gentle pulsing sound to remind you to ease off on the pressure.

## BrushSync replacement reminder



All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

## Ease into your new routine




Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

## The only true sonic technology




Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

## Up to 99.9% plaque removal



AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*\*

## Micro-droplet technology




Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

## High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

## Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

# Specifications

## Power

- Voltage: Multi-voltage charger

## Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

## Design and finishing

- Color: White with grey accents, White and Mint

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Battery indicator: Handle LEDs show charge level
- Handle: Slim ergonomic design
- Nozzle attachment: Easily snaps on and off
- Battery Life: 1-3 weeks\*\*\*\*\*

## Items included

- Charger: 1
- Brush heads: 1 G2 Optimal Gum Care

- Handle: 1 ProtectiveClean
- Nozzle: 1
- AirFloss Pro/Ultra handle: 1

## Cleaning performance

- Performance: Removes up to 99.9% plaque\*\*
- Health benefits: Cleans entire mouth in 60 sec, Improve gum health in 2 weeks\*
- Clinically proven results: As effective as floss, for gum health improvement, for plaque reduction, for reducing gum bleeding
- Customizable settings: Auto-burst or manual burst, Single/double/triple burst\*\*\*
- For best results: Change nozzle every 6 months

## Modes

- Clean: For exceptional everyday clean
- Gum Care: Gently massages gums

## Smart sensor technology

- BrushSync Replacement Reminder: Always know when to, replace brush heads



Issue date 2024-11-25

Version: 3.3.2

UPC: 0 75020 08384 5

© 2024 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* or your money back

\*\* when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.

\*\*\* From the treated areas; In a lab study, actual in-mouth results may vary

\*\*\*\* AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.

\*\*\*\*\* depending on burst setting used

\*\*\*\*\* than a manual toothbrush