

PHILIPS

AVENT

Pump, store, feed and care all-in-one set

Express, store, feed and care

Includes soft massage cushion

2 Bottles and 10 storage cups

Nipple protectors and cream

SCD292/01



More comfort, more milk* Pump, store, feed and care Single electric breast pump all-in-one set

The Breastfeeding Starter set is ideal for mums who plan to express their breast milk. It has all you need to pump and store your breast milk, and feed your baby. It includes speciality products to help you care for your breasts and nipples.

Express

- More comfortable pumping position due to unique design
- Features a gentle stimulation mode and 3 pumping settings
- Soft massage cushion with warm feel

Feed

- Unique Petals for a soft, flexible teat without collapse
- Natural latch on due to the wide breast-shaped teat

Store

- For use in fridge/freezer
- Easily pump in milk storage container

Care

- Help protect sore nipples during breastfeeding
- Silky-soft feel

Highlights

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means that you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow in the way that is most comfortable for you.

Soft massage cushion




Our massage cushion has a new soft, velvety texture that gives the skin a warm feeling for comfortable, gentle stimulation of your milk flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let-down.

Natural latch on



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the teat increase softness and flexibility without teat collapse. Your baby will enjoy a more comfortable and contented feed.

For use in fridge/freezer

The Philips Avent cups can be stored in the fridge or freezer and are dishwasher-safe

Easily pump in milk storage cup

With the milk storage cup adapter, you can easily pump in the milk storage container. No leakage, no spillage of your precious breast milk.

Help protect sore nipples during breastfeeding

Philips Avent Nipple Protectors are only designed for use when you have sore or cracked nipples and should be used with health professional advice.

Silky-soft feel



Silky-soft top-sheet and breathable, natural materials. Dermatologically tested.

Specifications

Design

Bottle design: Wide neck, Ergonomic shape

Breast pump design: Compact design

Country of origin

England

Material

Bottle: Polypropylene, BPA-free*

Breast pump: BPA-free*, Polypropylene

Teat: Silicone, BPA-free*

Breast pads: Dermatologically tested, Silky soft top-sheet, Natural breathable material

Storage cups: Polypropylene, BPA-free*

What is included

Milk storage cup lid: 10 pcs

Milk storage cup adapter: 2 pcs

Sealing disc for milk storage: 1 pcs

Natural bottle 125 ml/4 oz: 1 pcs

Extra soft Newborn Flow teat: 1 pcs

Natural bottle 9 oz: 1 pcs

Extra soft Slow Flow teat: 1 pcs

Breast pump body: 1 pcs

Standard-size cushion: 1 pcs

Base unit incl. tubing: 1 pcs

Milk storage cup (180 ml/6 oz): 10 pcs

Travel cover: 1 pcs

Breast pad sample packs: 1 pack (2 pcs)

Standard Nipple Protector: 2 pcs

Bra clip: 1 pcs

Moisturising Nipple Cream: 15 ml

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Settings: 1 Stimulation mode, 3 Expression settings

Breast care: Day, Night

Development stages

Stage: Pregnancy, 0-6 months

