

PHILIPS

sonicare

Sonic electric
toothbrush - Dispense

HealthyWhite

2 modes

2 brush heads

HX6782/34



Whiter, healthier teeth

Everybody loves a bright white smile. Let the natural whiteness of your patients' teeth shine through. HealthyWhite is proven to remove everyday stains in just two weeks, with the regular use of the Clean & White Mode.

Proven to improve oral health

- Naturally whiter teeth with patented sonic technology

Designed around you

- EasyStart

Guides you to follow dental recommendations

- Quadpacer
- Smartimer

Provides a superior clean

- Angled brush head neck for better reach to back teeth
- Safe on orthodontics, dental work and veneers
- Engineered to maximize sonic motion

Offers a customized brushing experience

- Clean and White Mode: proven to remove stains

Highlights

Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

Safe and gentle



The gentle motion of this Philips Sonicare power toothbrush is safe for dental work, including orthodontics, implants and veneers, so you can brush without worry.

Clean and White Mode



2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

Sonic technology




Philips Sonicare brush heads are vital to our core technology of high-frequency, high-amplitude brush movements performing over 31,000 brush strokes per minute. This dynamic cleaning action drives fluid deep between teeth and along the gumline, for a thorough yet gentle cleaning.

Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Quadpacer




Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

Smartimer



Two-minute timer helps ensure dental professional recommended brushing time.

EasyStart



Gently increases power over the first 14 uses to ease patients into the Sonicare brushing experience.

Specifications

Design and finishing

Color: Glacier Green

Power

Voltage: Multi-voltage charger

Technical specifications

Battery: Rechargeable

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

Brushing time: Up to 3 weeks**

Items included

Brush heads: 2 ProResults standard

Handles: 2 HealthyWhite

Charger: 1

Travel case: 1

Cleaning performance

Health benefits: Helps improve gum health

Speed: Up to 62000 brush movement/min

Performance: Better plaque removal than a manual toothbrush

Whitening benefits: Helps naturally remove stains

Timer: Quadpacer and SmarTimer

Modes

Clean: For exceptional everyday clean

White: Removes surface stains

