

Philips Zoom Chairside  
In office Kit

With Touch Up NW 22% CP-  
2pk

With Relief ACP Oral Care Gel

DIS569/01

## Philips Zoom In-Office Procedure Kits

For use with WhiteSpeed Lamp and Procedure gels

With the input of hundreds of Philips Zoom users, Philips has made administering the Philips Zoom whitening system easier, faster and safer than ever before.

**Impressive results fast**

- White bright smile to improve your patient's confidence

**Includes a Touch up treatment**

- Take home Nitewhite 22% CP

**Includes ACP Oral Care gel**

- Relief<sup>®</sup> ACP Oral Care Gel helps manage patient's sensitivity

# PHILIPS ZOOM!

# Specifications

## Consumer Trade Item

- Height: 8.64
- Width: 17.78
- Length: 31.75
- Net Weight: .8
- Gross Weight: .8
- GTIN: 08710103632047

- Country of origin: US
- Harmonized system code: 330690

## Outer Carton

- GTIN: 18710103632044

## Performance

- Whitening results: Maximum whiteness in 45 min.

# Highlights

## Includes ACP Oral Care gel

Recommendation is to place a light coating of Relief ACP (amorphous calcium phosphate) Oral Care Gel on the lingual surfaces of teeth on both arches before whitening treatment and to apply ACP Immediately following an in-chair procedure directly onto the teeth or in custom trays for 10 to 30 minutes

## Includes a Touch up treatment

For patients with severe tooth discolouration, a confident smile might not always come naturally. The innovative combination of Philips Zoom WhiteSpeed chairside whitening and Philips Zoom take-home treatment produces maximum results

## Impressive results fast

Your patients can get up to eight shades whiter teeth in 45 minutes, excluding preparation time



Issue date 2023-03-16

Version: 8.8.1

EAN: 87 10103 63204 7

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)