

# Where peace of mind meets efficiency

## Philips Emergency Care Device Readiness

Maintaining a fleet of monitor/defibrillators is a challenging but vital task. Caregivers rely on these devices for life-saving therapy – but since an emergency scenario could occur at any moment, each monitor/defibrillator must always be ready for use. This demands hospital biomedical engineers and EMS equipment managers stay fully informed on exactly which devices require maintenance at any given time.

To help meet these needs, we have developed the Device Readiness application. As part of the Philips Emergency Care Informatics Suite\* and specifically designed to work with the Philips HeartStart Intrepid monitor/defibrillator, this application is designed to give you the peace of mind that comes with knowing your Intrepid devices are ready for use when needed.

### Key benefits of Device Readiness


- Simplify device checks with integrated connectivity between Intrepid devices and Device Readiness
- Stay up to date on the readiness of your Intrepid fleet, with a customizable dashboard and configurable alerts
- Plan maintenance tasks and workloads effectively based on clear, accurate insights into device status
- Pave the way for a lower total cost of ownership for your Intrepid devices, by reducing the resources and effort needed to physically check each device

# An easy, efficient device readiness solution

## Discover the power of Device Readiness

Device Readiness is designed to help you maintain an overview of your fleet of Intrepid monitor/defibrillators. Data from Intrepid is automatically sent to a configurable dashboard which allows you to track key metrics. You receive up-to-date information on whether your equipment is ready for use without having to invest time and resources in manually inspecting each device for calibration and self-test data.

Ultimately, Device Readiness is designed to eliminate the hassle and inefficiency of maintaining a fleet of monitor/defibrillators - helping you to make sure that your devices are ready to help save lives.



## An integrated solution

Device Readiness is designed specifically to work with Intrepid monitor/defibrillators.

- Each Intrepid device can easily be configured to connect to your organization's Informatics Suite account
- Intrepid sends self-check data to the Informatics Suite cloud via WiFi or cellular
- Create a profile for your Intrepid, including name and assigned location plus up-to-date data transmitted from Intrepid

## Connectivity simplifies device checks

A customizable dashboard and configurable alerts keep you informed on the readiness of your Intrepid fleet.

- Hourly, daily and weekly self-test data feeds into your dashboard
- View shock button test status and last operational check date
- Receive alerts for NBP calibration
- Set alert preferences for maintenance or 'not ready for use' notifications

## Find and address devices that need maintenance


Self-check data helps reduce the resources needed to physically check each device.

- Review inspection note history
- Schedule inspections and EtCO<sub>2</sub> calibrations
- Access your dashboard via web browser
- Customize your dashboard and key performance indicators, filter information, and create an environment that supports your specific workflow

# When time and information are critical, our ALS solutions connect you... for life, where and when it matters.

To learn more about HeartStart Intrepid with Philips Emergency Care Informatics Suite, including Emergency Care Device Readiness, please contact your local Philips sales representative or visit our website at [www.philips.com/emergencycare-informatics-suite](http://www.philips.com/emergencycare-informatics-suite)

 **Emergency Care Informatics Suite**  
Philips Medical Systems  
22100 Bothell Everett Highway  
Bothell, WA 98021-8431, USA

 **HeartStart Intrepid**  
Philips Goldway (Shenzhen) Industrial Inc.  
Registered/Production address:  
No.2 Keji North 3rd Road,  
Nanshan District, Shenzhen, P.R. China 518057

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestraße 80, 20537  
Hamburg, Germany

EC REP  
CE 0123

