

PHILIPS

Emergency Care
Informatics Suite

AED Remote Monitoring



A smart service to ensure your AED is ready to use

AED Remote Monitoring

Automated external defibrillators (AEDs) have re-shaped the way we help protect our friends and colleagues against the risks of sudden cardiac arrests (SCAs). However, like any critical medical asset, AEDs require regular maintenance to help ensure they are ready to support you when you need them most.


Introducing Philips Emergency Care AED Remote Monitoring Suite, an all-in-one carefree ownership service

AED management challenges

There are many aspects to consider when managing AEDs:

- AED presence and accessibility
- Managing multiple AEDs at multiple sites
- Frequent inspections and preventive maintenance
- Compliance and mandatory updates
- Financial planning and device replacement management

How can you be sure your AED will be ready?





AED management redefined. It's as easy as 1...2...3!

Our AED Remote Monitoring Suite completely redefines the concept of managing AEDs. It is designed to be simple and easy to use with a carefree service subscription plan to help give you peace of mind. The AED Remote Monitoring Suite allows you to monitor the readiness of your AED assets, wherever they are. It utilizes a maintenance-free transmitter that connects your Philips AED to an automatic monitoring and notification portal.

1

Start with the industry standard, Philips HeartStart AEDs

Philips HeartStart AEDs

- Proven safe, reliable and effective¹
- Designed to be easy to use by trained or minimally trained users
- Fast, confident shock delivery < 8 sec²
- Daily checks, including pads functionality
- Ability to treat both children and adults in a variety of environments

2

Connect your AEDs, with the innovative SmartLink transmitter

Monitoring

- Remote monitoring; AED, pads and battery check
- Prompt notifications in case the AED needs attention
- Convenient status alerts
- Philips managed relationship with top cellular provider for cloud connectivity

3

Enjoy simple set-up, planning and ownership, with our remote monitoring service

Connected service

- Easy system for tracking parts that need replenishment
- Pre-paid subscription plan
- Downloadable snapshot of fleet readiness with one click
- Minimal visual inspections, enhancing sustainability

How AED Remote Monitoring could benefit you

Comprehensive support at every step, whatever your scale

Delivered on a subscription basis, our all-in-one carefree remote monitoring software provides you with a safe and reliable AED, an intuitive, subscription-based software platform and one-stop shop for device readiness tracking. Whatever your location or scale, we are by your side providing the confidence you need.

Key benefits

- Helps to be prepared if a sudden cardiac arrest (SCA) occurs
- Gain peace of mind and stay informed
- Focus on other priorities knowing your AEDs are being monitored
- High-quality, well maintained, 24/7 monitored AEDs
- Stay compliant with current local requirements
- Low total cost of ownership and predictable cost
- Improve sustainability by reducing unnecessary site inspections
- A truly carefree approach – we are your partner by your side



Are you ready to be
confident, connected and carefree?

Let's discuss how Philips can help you with your AED program.

www.philips.com/healthcare

Philips Medical Systems Nederland B.V.
Veenpluis 6
5684PC Best
The Netherlands



Philips Medical Systems
22100 Bothell Everett Highway
Bothell, WA 98021-8431, USA

1. Philips Medical Systems. (2009). Philips smart biphasic therapy. Retrieved from <https://www.usa.philips.com/healthcare/product/HC861304/heartstart-frx-automated-external-defibrillator>;
2. Nichol, G., Sayre, M. R., Guerra, F., & Poole, J. (2017). Defibrillation for Ventricular Fibrillation: A Shocking Update. *Journal American College of Cardiology*, 70(12), 1496-1509. doi:10.1016/j.jacc.2017.07.778;

© 2023 Koninklijke Philips N.V. All rights are reserved.
Philips reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.



00000086-00-00 * APR 2023
LC2272-025-009