



Removes up to 2X more plaque*

* vs a manual toothbrush

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

Proven to improve oral health

• Clinically proven safe and gentle

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

2 modes

- Clean: Standard mode for outstanding clean
- Massage

Offers a customized brushing experience

• Easy-start to acclimate to the Sonicare experience

Highlights

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

SmarTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Clear

Clean: Standard mode for superior teeth cleaning

Massage

- stimulates gums

Easy-start



Gently increases power over first 12 uses, to ease into the Sonicare experience



Specifications

Power

Voltage: 220-240 V, 50-60 Hz

Technical specifications

Battery: Rechargeable

Design and finishing

Color: Sky Blue

Service

Warranty: 2-year limited warranty

Ease of use

Battery indicator: Light in power button blinks to

indicate when to recharge

Brush head system: Easy screw-on brush heads

Brushing time: 20 2-minute brushings

Items included

Charger: 2 chargers Brush head holder Hard travel case: 2 **Cleaning performance**

Brush heads: 3 e Series std brush heads Performance: Better plaque removal than a

manual toothbrush

Health benefits: Helps improve gum health Speed: Up to 31000 brush movement/min

Timer: Quadpacer and SmarTimer

Whitening benefits: Helps remove tea, coffee

and other stains

© 2024 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2024-08-27 Version: 5.5.2

UPC: 0 75020 02710 8



