

Philips Sonicare Cordless  
Power Flosser 3000  
Oral Irrigator

**Quad Stream technology**

2 high performance nozzles: Quad  
Stream nozzle and Standard  
nozzle

Pulse Wave technology  
2 flossing modes and 3 intensity  
settings



HX3806/31

## Flossing reinvented with Quad Stream technology

Removes up to 99.9% of plaque\* in treated areas

Effortlessly thorough cleaning between teeth to improve your gum health. Unique X-shaped Quad Stream tip cleans more area with less effort for more effective flossing every time, so it's easy to get flossing right.

**Easy and convenient**

- Easy-to-fill 250 ml reservoir
- Effortless technique
- Standard tip for spot cleaning

**Faster, more effective flossing**

- 2 flossing modes, 3 intensities
- Complete clean in 60 seconds
- Pulse Wave technology guides from tooth to tooth
- Quad Stream technology for faster, more effective flossing

**PHILIPS**  
**sonicare**

## Oral Irrigator

Quad Stream technology 2 high performance nozzles: Quad Stream nozzle and Standard nozzle, Pulse Wave technology, 2 flossing modes and 3 intensity settings


# Highlights

### Quad Stream technology




Unique X-shaped Quad Stream tip separates flow into 4 water jets, which cover more area between teeth and along the gum line. Get a faster, deeper clean without the hassle of string floss.

### Standard tip for spot cleaning



A standard single-stream tip boosts pressure for for spot cleaning and removing food debris when you need it.

### Pulse Wave technology



Gentle pulses of water keep you going and guide you from tooth to tooth in Deep Clean mode — so you get it right every time.

### Complete clean in 60 seconds



Access those hard-to-reach areas with a nozzle that rotates 360 degrees, so that you can use the device in any orientation. Get a thorough 360° clean in just 60 seconds from start to finish.

### 2 modes, 3 intensities



Choose the clean that suits your smile. Clean mode delivers continuous flow for effective

HX3806/31


everyday cleaning. Deep Clean pulses for more thorough cleaning. Adjustable intensity level adapts for comfort.

### 250ml reservoir



Reservoir capacity of 250 ml holds enough water for a recommended 60-second clean. Simply twist and remove or use the side-fill door for easy refilling.

### Effortless technique



Just point, press and clean for more effective flossing. Aim the tip between your teeth and along the gumline — and let it do all the work for you!

## Oral Irrigator

Quad Stream technology 2 high performance nozzles: Quad Stream nozzle and Standard nozzle, Pulse Wave technology, 2 flossing modes and 3 intensity settings

# Specifications

### Power

- Voltage: Multi-voltage charger

### Service

- Warranty: 2-year limited warranty

### Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

### Performance

- Cleaning: Clean entire mouth in 60-90sec

### Items included

- Cordless Power Flosser: 1

- F1 Standard nozzle: 1
- F3 Quad Stream nozzle: 1
- Travel bag: 1
- USB charging cable: 1
- USB wall adaptor: 1

### Modes

- Clean: For exceptional everyday clean
- Deep Clean+: For an invigorating deep clean
- Intensities: 3

### Design and finishing

- Color: White



Issue date 2023-08-30

Version: 14.14.1

EAN: 87 10103 99180 9

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* In an in-vitro study. Actual results may vary.