

# PHILIPS

## sonicare

Sonic electric  
toothbrush

HealthyWhite+

2 modes

6 customizable settings

1 brush head

Whitening pen



HX8911/35



Removes up to 100% more stains in only 1 week\*  
\*than a manual toothbrush

Safely enhance your smile every day without giving up the staining foods and drinks you love! The Philips Sonicare Whitening rechargeable toothbrush is proven to remove up to 100% more stains for whiter teeth in only 1 week\*.

#### Proven to improve oral health

- Improves gum health in only two weeks
- Removes up to 7x more plaque than a manual toothbrush
- Removes up to 100% more stains than a manual toothbrush\*

#### Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

#### Provides a superior clean


- Philips Sonicare's advanced sonic technology
- DiamondClean brush head for your best whitening

#### Offers a customized brushing experience

- 6 customizable settings
- Clean and White Mode: proven to remove stains

# Highlights

## Whitens teeth




This Philips Sonicare toothbrush removes up to 100% more stains for whiter teeth in just 1 week.

## DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles that are medium stiff to effectively yet gently scrape away plaque. This Philips Sonicare electric toothbrush provides an exceptional cleaning and whiter teeth compared to a manual toothbrush.

## Clean and White mode



2 minute Clean mode offers the ultimate in plaque removal (Dental professional recommended brushing time). White mode removes surface stains, brightens and polishes

teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.

## 6 customizable settings




Philips Sonicare Stain Removal toothbrush allows you to mix and match modes (Clean mode, White mode) and 3 intensities for maximum comfort and your best clean.

## Removes up to 7x more plaque



Removes up to 7 times more plaque from hard to reach places than a manual toothbrush.

## Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks.

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

## SmarTimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## The only true sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.



# Specifications

## Design and finishing

**Color:** Frost White

## Power

**Voltage:** 110–220 V

## Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** Up to 3 weeks\*\*

**Battery type:** Lithium ION

## Service

**Warranty:** 2-year limited warranty

## Ease of use

**Battery indicator:** Light shows battery status

**Brush head system:** Easy click-on brush heads

**Handle:** Slim ergonomic design

**Brushing time:** Up to 3 weeks\*\*

## Items included

**Handles:** 1 HealthyWhite+

**Charger:** 1

**Travel case:** 1

**Brush heads:** 1 DiamondClean standard

**Whitening toothpaste:** 1

**Whitening pen:** 1

## Cleaning performance

**Performance:** Removes up to 100% more plaque

**Health benefits:** Improves gum health in only two weeks

**Speed:** Up to 62000 brush movement/min

**Timer:** Quadpacer and SmartTimer

**Whitening benefits:** Whitens up to 2 shades in 2wks

## Modes

**3 intensities:** Low, Medium, High

**Clean:** For exceptional everyday clean

**White:** Removes surface stains

