


Philips Sonicare
Electric Toothbrush &
Dental Flosser

w/ ProtectiveClean Electric
Toothbrush
AirFloss Pro
1 brush head & 1 nozzle



HX8491/74

Healthier gums in 2 weeks, guaranteed*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra**** is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is clinically proven to be as effective as floss for gum health.**

Improves oral health

- Clinically proven to be as effective as floss for gum health**
- Up to 99.9% plaque removal***
- Helps prevent cavities between teeth

Optimise your brushing

- Choose from three modes

Superior technology

- Air and micro-droplet technology
- High-performance nozzle
- Triple-burst customisable settings
- Connects smart brush handle and smart brush heads


The easiest way to effectively clean between teeth

- Takes just 60 seconds: Point, Press, Clean!
- An easy way to start a healthy habit
- Total daily confidence with a fresh and clean mouth

PHILIPS
sonicare

Highlights

Three modes




With a choice of three modes, this toothbrush allows you to customise brushing to your own needs. Clean mode is standard for superior cleaning. White is the ideal mode for removing surface stains. And Gum Care mode adds an extra minute of reduced-power brushing, so you can gently massage your gums.

BrushSync technology




A microchip-enabled technology that detects and synchronises the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.**

Micro-droplet technology




Our clinically proven results are possible thanks to our unique technology, which combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss Pro/Ultra is an easy way to clean deeper between teeth, helping to form a healthy habit.

Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple) and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Fresh and clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, and then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

Specifications

Power

- Voltage: Multi-voltage charger

Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

Design and finishing

- Colour: Black

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle: Slim ergonomic design
- Battery indicator: Handle LEDs show charge level
- Nozzle attachment: Easily snaps on and off
- Battery Life: 1–3 weeks****

Items included

- Charger: 1
- Brush heads: 1 W2 Optimal White
- Handle: 1 ProtectiveClean
- AirFloss Pro/Ultra handle: 1
- AirFloss Pro/Ultra nozzle: 1

Cleaning performance

- Health benefits: Improve gum health in 2 weeks*, Cleans entire mouth in 60 secs
- Performance: Removes up to 99.9% plaque**
- Clinically proven results: As effective as floss, for gum-health improvement, for reducing gum bleeding, for plaque reduction
- Customisable settings: Auto-burst or manual burst, Single/double/triple burst***
- For best results: Change nozzle every 6 months



Issue date 2023-05-03

Version: 5.0.2

EAN: 87 10103 87901 5

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* or your money back

** when used in conjunction with a manual toothbrush and anti-microbial mouthwash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under the Support tab for further details.

*** From the treated areas; In a lab study, actual in-mouth results may vary

**** AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.

***** depending on burst setting used