

# PHILIPS

Auto CPAP  
Advanced

DreamStation 2

Includes Ramp Plus

Color touchscreen

Integrated humidifier

3 adjustable comfort settings

HH1551/00



## Advanced care designed to meet your needs

### Includes 15mm Standard Tubing

DreamStation 2 Auto CPAP Advanced is designed with your comfort and ease in mind. Now with Ramp Plus, you can pick a starting pressure that is comfortable for falling asleep. It also features a colorful touchscreen for easy navigation.

#### **DreamStation 2**

- Designed for comfort
- Easy to use
- Quiet sound
- Small size
- Keeps you motivated

# Highlights

## Comfort



With Ramp Plus, you have the ability to pick the starting pressure that is comfortable to fall asleep to and will adjust should an event occur.

## Easy



Three comfort settings are designed to be easy to adjust while therapy is running. Integrated connectivity within the device allows you to stay connected with your care team. 93% of users describe the CPAP as “easy to use”.<sup>1</sup>

## Quiet



Users rated DreamStation 2 Auto CPAP Advanced quieter than their prescribed device.<sup>2</sup>

## Small - 15mm tubing



DreamStation 2 Auto CPAP Advanced is 22% smaller and 16% lighter than the other leading CPAP brands.<sup>3</sup>

## DreamMapper



DreamMapper mobile app allows you to track your progress, set goals and get helpful tips. People who struggled with sleep apnea therapy had a 283% higher success rate when they used DreamMapper.<sup>4</sup>

# Specifications

## Product details

**Weight:** 2.3 lbs without power supply; 2.94 lbs with power supply

**Dimensions:** 10.11 in L x 6.24 in W x 3.34 in H

**Package contents:** Auto CPAP Device

**Package contents:** 1 Disposable Ultra-fine filter, 1 Reusable pollen filter, 15mm Tubing, Carry bag, Power cord, Power supply, User Manual

**Tubing compatibility:** 12mm standard/heated, 15mm standard/heated, 22mm standard

## Technical Specifications

**Pressure range:** 4 to 20cm H<sub>2</sub>O

**Flex Pressure Relief:** Off, 1, 2, 3

**Ramp time:** 15 minutes, 30 minutes, 45 minutes

**Humidification:** Adaptive humidification and heated tube settings (optional)

**Altitude compensation:** Automatic

**Filters:** Reusable pollen; disposable ultra-fine

**Device controls:** LCD color touchscreen; push therapy button

**Electrical requirements:** 100–240 VAC, 50/60 Hz, 1.0–2.0 A

**Modes:** CPAP, Auto CPAP, Fixed CPAP



\* 1. 30-day 2020 patient preference study with N=444 users. Prescribed devices included Philips DreamStation and ResMed AirSense 10.

\* 2. Based on users' rating of sound and comparison from 30-day 2020 patient preference study with N=444. Prescribed devices included Philips DreamStation and ResMed AirSense 10.

\* 3. DreamStation 2 Auto CPAP Advanced is 22% smaller and 16.4% lighter than ResMed AirSense 10.

\* 4. Based on 2020 patient preference study with N=444.