

# PHILIPS

## AVENT

Manual breast  
pump

with VIA Storage System

SCF300/12

 Product Picture

## Express more. Store easily.

Express, store and feed in easy steps

The unique Philips Avent Manual Breast Pump SCF300/12 offers you total control of the suction level and pumping rhythm. It assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump.

**Benefits and features:**

- Gently stimulates fast milk flow

**Clinically proven effective**

- Patented Let-down Massage Cushion

**Fits all Philips Avent breast pumps and teats**

- Use adapter to express and feed

**For safe storage and transport**

- Leak-proof twist-on system

## Highlights

### Gently stimulates fast milk flow

As you pump, the patented Philips Avent Let-down Massage Cushion flexes in and out, gently massaging the area around the nipple to help stimulate natural milk flow.

### Patented Let-down Massage Cushion



The five soft petal massagers of the Philips Avent breast pump gently flex in and out as you pump and work with the vacuum to imitate your baby's suckling.

### Leak-proof twist-on system

For safe storage and transport

### Adapter to express and feed

Use the adapter to express milk and feed with the same cup.

## Specifications

### Country of origin

England

### Development stages

Stage: 0–6 months

### Material

BPA-free\*

### What is included

Milk storage cup lid: 5 pcs

Milk storage cup adapter: 2 pcs

Extra soft Newborn Flow teat: 2 pcs

Newborn nipple travel pack: 2 pcs

Spare parts: 1 pcs

Pre-sterilised VIA Cup (180 ml/6 oz): 5 pcs

Manual breast pump: 1 pcs

