




Philips Sonicare AirFloss
Interdental - Rechargeable

w/ 3 Nozzles

1 Nozzle holder
Travel case



HX8255/02

Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is the easiest way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to remove plaque where brushing can't.

Improves oral health

- Removes more plaque between teeth
- Helps improve gum health in just two weeks
- Helps prevent cavities between teeth

Achieve a 30-second clean

- Cleans your entire mouth in just 30 seconds
- Point. Press. Clean!

Easy to use

- Air and micro-droplet technology
- Comfortable ergonomic handle
- Guidance tip ensures correct placement


Helps build a healthy habit

- An easy way to improve cleaning between teeth

PHILIPS
sonicare


Highlights

An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

Removes more plaque between teeth



Removes more plaque between teeth

30-Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from

space to space. It will automatically burst every second.

Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

Air and micro-droplet tech



AirFloss utilises unique air and micro-droplet technology to send a burst of pressurised air and water or mouthwash between your teeth. The one-touch activation button is as simple as a mouse click.

Ergonomic handle

Designed to be easy to manoeuvre in tight spaces, AirFloss has a comfortable ergonomic handle and a non-slip grip.

Guidance tip



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

Specifications

Power

- Voltage: Multi-voltage charger

Technical specifications

- Operating time: Replace nozzle every 6 months
- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

Design specifications

- Colour(s): White with green accents

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle: Slim ergonomic design

- Battery indicator: Flashing indicates low battery
- Battery Life: Up to 2 weeks
- Nozzle attachment: Easily snaps on and off

Items included

- Travel case: Soft travel case
- Charger: 1 multi-voltage charger
- AirFloss Handle: 1
- AirFloss Nozzle: 3
- Nozzle holder: 1

Cleaning performance

- Auto-burst feature: Cleans entire mouth in 30 secs
- Health benefits: Improves gum health in only two weeks
- Nozzles: Best results change every 6m
- Performance: Removes up to 5x more plaque*



Issue date 2023-05-19

Version: 5.0.2

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* than manual tooth-brushing alone