

PHILIPS

sonicare

Sonic electric
toothbrush

DailyClean

1 mode

1 brush head



HX3411/04



Clinically proven to provide superior cleaning ^{*}

Up to 3x more plaque removal*

Designed to be the perfect toothbrush to transition from your manual toothbrush. Until now, it feels like you've never brushed your teeth before.

Proven to improve oral health

- Removes up to 3x more plaque than regular brushing

Designed around you

- Easy start program builds your Philips Sonicare routine
- Sonicare rechargeable toothbrush
- 2 minute brushing time recommended by dentists

Innovative technology

- Sonic advanced technology removes food and plaque

Optimize your brushing

- Battery life
- Click-on brush head system
- Just one click to change brush heads

Highlights

Sonicare 1100 — a deeper clean



For an exceptional daily clean, Sonic technology pulses water between your teeth and gums to break up plaque and blast away food remnants.

Smart functionalities




With a gentle pulse reminds you to change to a new part of your mouth, and signals when 2 minutes of brushing time is up. This sonic electric toothbrush gives you the most optimal brush, every time.

Ease into your new routine



Starting something new can take some getting used to. Our easy-start program gives you the option of a gradual, gentle increase in brushing power with the first 14 uses with your new toothbrush.

Battery life



Lasts up to 14 days between charges.

Click-on brush head system



This brush head clicks on and off your brush handle for a secure fit and easy maintenance and cleaning. It fits Philips Sonicare toothbrush PowerUp Battery, Essence+/Elite+/CleanCare+ handles.

Say goodbye to plaque



The Philips Sonicare 1100 gets at all those hard-to-reach areas that a manual toothbrush can't get to, with contoured bristles designed to make brushing easy.

Use at home or on-the-go

The long-lasting battery means it lasts up to 14 days between charges. Perfect for on the go!

Easy to maintain



Dental professionals recommend changing your brush head every 3 months. This brush head clicks on and off your brush handle for a secure fit and easy maintenance and cleaning.

Specifications

Power

Voltage: 110–220 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 14 days

Battery type: NiMH

Design and finishing

Color: Mint

Service

Warranty: 2-year limited warranty

Ease of use

Brush head system: Easy click-on brush heads

Handle: Ergonomic design

Items included

Charger: 1

Brush heads: 1 C1 SimplyClean standard

Handle: 1 DailyClean

Cleaning performance

Performance: Removes up to 3x more plaque*

Speed: Up to 31000 brush movement/min

Timer: Quadpacer and SmarTimer

Modes

Clean: For exceptional everyday clean

